

Aesculap[®] activ C

Cervical Disc Prosthesis
Retain Mobility



Aesculap Spine

activ C

Back to life enjoyed to the full



- † Natural mobility
- † During surgery and in everyday life
- † Stability
- † The step into an eventful future supported by Clinical evidence

activ C

V A T I O N



Natural mobility



+ Retaining mobility of your patients

I T Y



Physiological center of rotation



■ COR in central position

■ Close to COR of degenerated discs according to latest study results

Suchomel et al.: Sagittal segmental alignment after activ C cervical arthroplasty after 1 year follow-up, Abstract #35, ISASS 2012, Barcelona.

Penning L.: Normal movements of the cervical spine. Am J Roentgenol. 1978;130:317-26.

Dvorak J. et al.: In vivo flexion/extension of the normal cervical spine J of Orthopaedic Research, 1991;9(6):828-34.

Reconstruction of lordosis



■ Good balancing and alignment

■ Significant correction of segmental angulation

Suchomel et al.: Does sagittal position of the CTDR related center of rotation influence functional outcome? Prospective 2 year follow-up analysis, Abstract #35, ISASS 2012, Barcelona.

Sustainable restoration of height



■ Low profile

■ Good contact to bone surface

■ Anatomical adapted footprint and shape

■ No subsidence

■ No dislocation


■ Restoration of Height

Suchomel et al.: Does sagittal position of the CTDR related center of rotation influence functional outcome? Prospective 2 year follow-up analysis, Abstract #35, ISASS 2012, Barcelona.

Meisel et al.: Does CTDR have a lower risk of device subsidence compared to ACDF? 2 year results of a prospective multi-center study, Abstract #249, ISASS 2012, Barcelona.

activ C

During surgery and in everyday life



+ Multilevel treatment through intelligent implant design

T Y



Multilevel treatment



Combination of spikes and keel

Treatment of 2 or 3 levels is possible. Prevention of vertebral body split caused by double keel prosthesis

Safe midline positioning



Keel on the inferior side

Solid anchorage in the inferior vertebral body and accurate positioning in the midline

Safe preparation of the keel bed



Intelligent instrumentation and reamer guidance

No chiseling

Easy and gentle preparation of the keel bed; risk of spinal cord and blood vessel injuries reduced

activ C

Stability in every situation



STABI

+ High stability and accurate placement

L I T Y



Anatomic design



Large contact area through convex shape of the superior prosthesis plate

Guided by the anatomy of the vertebral body, allowing accurate placement and fitting
▶ primary stability

Solid anchorage cranial + caudal



Spikes on the superior plate

Improved anchorage in the dome - superior to keel anchorage ▶ primary stability

Keel on the inferior plate

Stable anchorage in the inferior vertebral body and solid support against rotation and lateral dislocation ▶ primary stability

Rapid and safe osteointegration



Plasmapore® coating

Fast ingrowth of bone cells
▶ secondary stability



Clinical evidence for an effective treatment and a better quality of life



Qualitative evaluation of study device



- ✓ All devices show ideal lateral placement (< 2 mm of midline placement in the m/l-direction)
- ✓ All devices intact (no device disassembled, loose or fractured)
- ✓ No device subsidence (≥ 3 mm) observed
- ✓ No device migration (> 3 mm) observed
- ✓ No device expulsion (≥ 50 % of the a/p dimension of the device extends beyond the anterior margin of the disc space)

Clinical Results - NDI & VAS





- Reduction of NDI over all time points
- Significant reduction of NDI from pre-OP to 6 weeks and from 6 weeks to 6 month
- Significant reduction of VAS Neck & Arm Pain from pre-OP to 6 weeks and from pre-OP to 1 year

Functional Results - ROM & Disc Height



- Slight correction of degenerative instability directly post-OP
- Very good maintenance of ROM over all time points
- Significant restoration of disc height and maintenance of disc height from 6 weeks post-OP (after spike penetration) to all other time-points

Functional Results - Correction of Disc Angle / Lordotization



- Significant restoration of disc angle with significant post-OP lordosis
- Maintenance of lordosis over all time points
- Significant correlation of COR of implant and the extent of achieved lordotization (the more anterior the COR-I, the better the correction of disc angle)

