




PARADIGM SPINE
the movement in spine care



Dynamic Stabilization System **DSS**™

Lumbar Dynamic Stabilization



DSS™

Dynamic St

The basic concept of the DSS system was developed in cooperation with the Institute of Orthopedic Research and Biomechanics of the University of Ulm.

The goal was the development of a dynamic connector that reduces the natural range of motion (ROM) by approximately 50% whilst maintaining the physiological center of rotation.

To maintain physiological kinematics length adaptation under flexion / extension movements is mandatory. Based on a validated Finite Element model of a lumbar spinal segment, the mechanics of the dynamic coupler was developed and confirmed in a biomechanical experiment. The results were published in SPINE by Prof. Wilke.*

The DSS system was launched in 2008 after 2 years of development.

The product obtained FDA approval in 2009.

In contrast to the 1st generation dynamic systems, the DSS system has following features:

- Controlled movement of the dynamic coupler
- Combination of fusion and dynamic stabilization
- Controlled dynamic stabilization of hypermobile segments while effectively protecting adjacent levels above and below fusion
- Restoration of physiological center of rotation
- Shock absorption; axial compression
- Control of rotation and stress reduction on facet joints
- Significant reduction of shear forces

* Prof.H-J. Wilke, Frank Heuer, PhD, Hendrik Schmidt, PhD
2009, Prospective Design Delineation and Subsequent In Vitro Evaluation of a New Posterior Dynamic Stabilization System, SPINE Volume 34, Number 3, pp 255–261, Institut für Orthopädische Forschung und Biomechanik, Ulm, Deutschland

Next generation dynamic stabilization

Progressive degeneration of a spinal motion segment leads to an increase of the range of motion in flexion / extension, lateral bending and rotation. The neutral zone is increased as well. Hypermobility causes back pain and muscle spasms.


The DSS system achieves effective stabilization by reducing the segmental mobility (by the stiffness of the coupler), as well as significant reduction of rotational and translational movements (protection of the facet joints).

Extreme segmental movements are eliminated. The range of motion is limited to an "active zone".


In cooperation with different anatomical institutes Paradigm Spines offers DSS training courses focusing on the paraspinal approach. Clinical visits in various reference centers are also possible.

For scientific documentation an internet-based registry is available, that allows for easy documentation of patient data and allows easy evaluation and exportation of charts and graphics for PowerPoint presentations, patient conversations or publications.


abilization




Reduction of Range of Motion ¹



Maintenance of Center of Rotation ²




Pedicle Displacement ²




¹ Prof.H-J. Wilke, Frank Heuer, PhD, Hendrik Schmidt, PhD
2009, Prospective Design Delineation and Subsequent In Vitro Evaluation of a New Posterior Dynamic Stabilization System, SPINE Volume 34, Number 3, pp 255–261, Institut für Orthopädische Forschung und Biomechanik, Ulm, Deutschland

² Bryan W. Cunningham, MSc, Dennis Collieran, BS, Gwen E. Holsapple, BS, Karen A. Adams, BS, Paul McAfee, MD (2006): Lumbar Spine Kinematics – a Radiographic Assessment of the Change in Interpedicular Distance throughout the Range of Motion. SAS 6 Global Symposium on Motion Preservation Technology, Montreal. Orthopaedic Spinal Research Laboratory, St. Joseph Medical Center, Towson, Maryland, USA.

Modular Design



Modularity of the system allows combination of fusion and dynamic (hybrid) solutions. Polyaxial positioning of $\pm 14^\circ$ possible through the coupler-spacer interface.



Biomechanic

Stabilization

Pain Relief

- Controlled dynamic stabilization of hypermobile segments
- Restoration of physiological center of rotation
- Shock absorption; axial compression
- Control of rotation and stress reduction on facet joints

Hybrid Solutions

- Modularity of system allows combination of fusion and dynamic stabilization
- Adjacent levels are effectively protected

Ease of Use


- Cannulated screws support minimally invasive application
- Paraspinal muscle splitting technique minimizes iatrogenic trauma
- Easy to use instrumentation

DSS™ Implant – Indications

- The DSS system is indicated for skeletally mature patients with degenerative disc disease (DDD) at one to three levels from L1 to S1, including conditions up to grade 1 spondylolisthesis.
- DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies, e.g.:
 - Segmental sagittal hypermobility
 - Position of segment in X-Ray (Standing Position)
 - Disc height (Black Disc)
 - Modic signs (MRT)
 - Painful degeneration of facet (positive infiltration)
 - Painfree response from segmental epidural injection

Modular





Semi-rigid stabilization + Fusion

Fusion



System

DSS™


Modular System

- Modularity of system allows combination of fusion and dynamic stabilization.
- Pathologies can be addressed individually by motion segment.
- Adjacent segments are protected.

Surgical Tech


1. Patient Positioning

- The patient is placed in neutral position.



2. Paraspinal Intermuscular Approach

- Intermuscular approach between multifidus and longissimus muscle (along the fascia).
- Cannulated screws in combination with K-wires support minimally invasive, tissue-sparing application.



3. Pedicle Preparation

- After standard skin incision via a midline or bilateral paramedial approach the pedicles are exposed and the pedicle entry point is marked under X-ray control using the pedicle awl.
- This awl opens the cortical bone and creates an initial hole in the pedicle as a starting point for further pedicle exploration and drilling.



Technique 1

DSS™


4. Screw Insertion

- Cannulated, ratcheting T-Handle connects to self-capturing Pedicle Screw Inserter.
- Facilitates rapid and safe screw insertion.
- Screw held securely - control during insertion.
- Cannulated screws in combination with K-wires support minimally invasive, tissue-sparing application.



5. Pedicle Screw Implant Guide

- Securely fixed to screws facilitate implant delivery for true MIS benefits.
- The Extension Rod (DAT40410) has a predetermined breaking point.



6. Delivery of Spacers and Washers over the Implant Guides

- The spacers and washers are delivered over the implant guides.
- The pusher slides over the implant guides.



Surgical Tech


7. Caliper Adapters

- Depending on the desired coupler position (high or low assembly) the caliper adapters are attached over the extension rods.



8. Measuring

- The caliper measures the distance between the caliper adapters to determine the exact coupler length.
- Caliper adapters face one another.
- Caliper is positioned between adapters and allowed to expand.
- Exact coupler length is read on graduated scale.



9. Coupler Selection

- Choose appropriate coupler according to caliper measurements.
- Adjust coupler length using the Coupler Length Adjuster.
- Note: one full turn represents 0.5mm length change.



10. Coupler Positioning

- The coupler pusher seats the coupler onto the construct.
- In the low assembly the spacer is used to maintain the coupler in place.
- Compression and distraction Instruments may be utilized to bring the slotted Coupler down.



Technique 2

DSS™

11. Nut Placement

- Slide breaking nut over the extension rod (shorter side down) and pre-tighten nuts using the nut driver.
- During nut tightening the retaining forceps slotted coupler or the retaining forceps rigid coupler should be utilized.
- Pre-tighten nuts using the nut driver.



12. Breaking Nuts and Final Tightening

- Ratchet goes on top of the nut driver.
- The upper part of the nut breaks off at a safety torque of 12 Nm.
- Final nut tightening using the nut driver and the counter torque instrument.
- Care has to be taken to keep the dynamic coupler aligned during tightening.



13. Final Position




Assembly

DSS™

Monosegmental



Bisegmental




Instructions


Three Level



Washers („blue“) are placed at the end points of assembly.
Note: If a hybrid construct is assembled, the fusion coupler has to be placed first.



Washers („blue“) are placed at the end points of assembly.
Note: If a hybrid construct is assembled, the fusion coupler has to be placed first.




Patient Cases

DSS™



Preoperative



Case 1: Male, 38 years

- Low back pain for 8 years with occasional leg pain while walking. ODI 64%, VAS back pain 7 and leg pain 3.
- Diagnosis: Chronic lumbago with DDD and medial protrusion.
- Previous therapy: Failed conservative treatment including facet infiltration.
- Surgery: Dynamic stabilization with DSS™ at level L4/5.
- Follow-up after 6 months: Reduction of ODI to 25% and improvement of VAS back pain 3 and leg pain 1.





Postoperative






Preoperative



Postoperative





**Case 2:
Male, 52 years**

- Low back pain for 6 years. ODI 63%, VAS back pain 7 and leg pain 3.
- Diagnosis: Severe DDD with rotational instability at level of L4/5, osteophytes left as signs for rotational instability and initial degeneration in L3/4 with hypermobility.
- Previous therapy: Failed conservative therapy for 5 years with constant pain medication. Since two years physiotherapy without improvement of symptoms.
- Surgery: Dynamic stabilization with DSS™ at level L3/4 and fusion at level L4/5 with additional implantation of TLIF at L4/5.
- Follow-up after 6 weeks out: Reduction of ODI to 34% and improvement of VAS back pain 4 and leg pain 1.

Patient Cases

DSS™


Preoperative



Case 3: Female, 52 years




- Low back pain for many years during standing and walking, virtually no physical activity possible through back pain. ODI 66%, VAS back pain 7 and leg pain 4.
- Diagnosis: Isthmic spondylolisthesis L5/S1, facet arthrosis L4/5 grade 3 (Fujiwara) and functional retrolisthesis L4/5.
- Previous therapy: Failed conservative treatment including continuous physiotherapy for 5 years. Constant pain medication according to WHO1 classification for 3 years.
- Surgery: Dynamic stabilization with DSS™ at level L4/5 and fusion at level L5/S1 with additional implantation of TLIF at L5/S1.
- Follow up: Full reduction and normalization on L4/5, L5/S1. Direct postop drop of back pain to VAS 5. After 6 months VAS back pain 3 and no leg pain.

Postoperative






Preoperative



**Case 4:
Female, 64 years**

- Low back pain for 10 years with occasional pseudoradicular leg pain. Walking distance 50m, ODI 68%, VAS back pain 8 and leg pain 5.
- Previous therapy: Failed conservative treatment, disc protrusion L3/4 to L5/S1 with decompression. Spondylarthrosis grade 3 as well as scoliosis
- Surgery: Dynamic stabilization with DSS™ at level L3/4 and fusion at level L4-S1 with additional implantation of TLIF at L4/L5 and L5/S1.
- Follow-up after 6 months: Improvement of back pain from VAS 8 to 4. Walking distance increased to 250m. ODI 38%.
- Follow-up after 18 months: ODI 6%, VAS back pain 2 and leg pain 2.


Postoperative



Objective

During the last several years, dynamic stabilization systems in spine surgery have become important clinical options for surgeons and patients. However, there are few long-term studies which present the safety and effectiveness of these new methods in comparison to the gold standard. Randomized clinical studies are often long, very complex and costly.

As an alternative to these studies, Paradigm Spine established a worldwide Registry to support prospective data collection on safety and efficacy of our DSS™ technology. The data will be documented in the Registry which is available exclusively through Paradigm Spine.



Description


The Paradigm Spine Registry is an internet-based system for patient documentation. Relevant clinical data will be collected and standard reports are possible.

The data entry is conducted online:

- Preoperative forms (clinical and radiographic evaluation)
- Surgery details (diagnosis and treatment)
- Postoperative forms (clinical evaluation, different follow-up intervals are possible)
- Patient forms (VAS and ODI resp. NDI scores, pre- and postoperative)
- Hospital Discharge details (postoperative treatment)

The Registry is available for documenting the DSS™ system and the DCI™ implant.

Paradigm
www.paradigm.com





Reports

The Registry provides several clinical reports. The reports are up to date at any time and can be easily generated via mouse click. Among others, the following reports are available:

- VAS score
- ODI resp. NDI index
- SF-12v2® score

Reports may be generated for a single patient, as well as for all patients of one site. In addition to that a comparison against the entire data of the Registry is possible. The graphs can be easily transferred into various IT programs (e.g. PowerPoint) and used for presentations.



Paradigm Spine Registry
digmspine.com/registry

Advantages

The Paradigm Spine Registry offers several benefits to its users:

- Patient documentation online
- Control of clinical results
- Internet-based database, less paperwork
- Platform for studies and user groups
- Fast clinical reports and statistics
- Eases preparation of lectures and podium presentations
- Basis for publications
- Worldwide access

Data Privacy

The data privacy regulations have been followed, and all data collection is anonymous. The identification of a patient is only traceable by the surgeon.


Access to the Registry is only available with a separate personal password.

Instruments

DSS™




Additional instruments



① Pedicle preparation and screw insertion



② Measuring




③ Implant positioning / tightening





Pedicle Screw Sets

(includes 2 cannulated Pedicle Screws per package)

	Length	Ø 6	Ø 7	Ø 8
	35 mm		DAI47035	DAI48035
	40 mm	DAI46040	DAI47040	
	45 mm	DAI46045	DAI47045	DAI48045
	50 mm	DAI46050	DAI47050	
	55 mm	DAI46055	DAI47055	DAI48055



Closure Sets

(includes 2 Nuts, 2 Spacers, 2 Washers per package)



	Standard (Breakoff)	Optional (Non-Breakoff)	
	DAI06041	DAI06042	

Coupler

(1 coupler per package)

	Dynamic	
	Sizes	Ref. / Order No.
	24-26 mm	DAI35024
	26-28 mm	DAI35026
	28-31 mm	DAI35028
	31-35 mm	DAI35031
	Fusion	
	Sizes	Ref. / Order No.
	20 mm	DAI70020
	22-25 mm	DAI70022
	25-31 mm	DAI70025
	31-43 mm	DAI70031

Optional

	Spacer with Thread
	DAI06020
	Implant Guides, Sterile Packed
	DAT40611



Biocompatibility

All DSS™ implants consist of titanium 6-aluminium 4-Vanadium (ISO 5832-3).

MRT

Titanium is a non ferromagnetic material, therefore magnetic resonance imaging can be done.

Retraceability

All DSS™ Implants are delivered sterile packed.


Patent

Patent applied for DSS™ system.


Implants

DSS™


Pedicle Screw Set



Dynamic Coupler



Fusion Coupler





PARADIGM SPINE

the movement in spine care

Paradigm Spine GmbH
Eisenbahnstrasse 84
D-78573 Wurmlingen, Germany

Tel +49 (0) 7461 - 96 35 99 - 0
Fax +49 (0) 7461 - 96 35 99 - 20

info@paradigmspine.de
www.paradigmspine.com

Product not available in the USA