





GLOBAL[®]
biomedica



The Lumbar BMD
Titanium PLIF GW Implant[®]

Surgical Technique Guide



3D-printed spinal fusion Global BMD Titanium Implants®

TABLE OF CONTENTS

Preparation and Approach	4
Decompression and Discectomy	4
Distraction	5
Determination of Implant Size	5
Implant Insertion	6
Completion of Surgery	6
Implant Removal	7
List of Implants	8
Tool List	11
Important Information on the Global BMD Titanium PLIF GW Implant®	12

The following general Surgical Technique Guide is for illustrative purposes only. As with all surgical procedures, the technique used in each case will depend on the surgeon's medical judgment as to the best treatment for each patient. Only those individuals with specialized training and experience in spinal surgery should attempt to use the BMD Titanium PLIF GW Implant®. Detailed preoperative clinical and diagnostic evaluation followed by carefully executed surgical technique is essential. Refer to the Instructions for Use (IFU) for a complete list of prescribing information. This technique guide was developed in conjunction with health care professionals. Additional information about the medical device is available from Global Biomedica s.r.o.



Figure 1
Patient positioning

Preparation and Approach

The patient is placed on the OR table using the standard positioning in cases of Posterior Lumbar Interbody Fusion (Figure 1). Make sure the abdomen is free positioned to avoid pressure on the large vessels and to minimize blood loss. A radiolucent OR table is recommended, as X-ray shall be used to confirm identification of the affected disc and in a later stage the position of the implant. Mark the affected segment after fluoroscopy control. A standard incision is performed over the level(s) to be instrumented. Expose the spinous process and facet joint on the affected side. Use a high-speed drill, or Kerrison rongeur Resect the necessary part of lamina and facet joint to approach the spinal canal, avoid the completely resection of facet joint, if possible.



Figure 2
Endplate preparation

Decompression and Discectomy

To ease the discectomy procedure, insertion of pedicle screws and performing a slight distraction over the pedicle screws prior to the discectomy can be advantageous. The neural foramen and spinal canal are decompressed as necessary. Dural sac is gently retracted medially using a nerve root retractor. The posterolateral portion of the annulus is then exposed and a window is created to gain access to the intervertebral space. The dural sac and neural structures are continually protected by the nerve root retractor. After the disc is cleared, the endplate cartilage should be removed carefully, leaving the upper and lower bony endplate intact (Figure 2). Injury to the bony endplates may lead to implant subsidence.



Figure 3
Intervertebral disc space

Distraction

An adequate distraction of the intervertebral disc space is one of the preconditions for the primary stability required after cage implantation. For this purpose Distractors are available in several heights (Figure 3). Each size of Distractor has grooves, that indicate the length of 20, 24, 28 and 33 mm respectively. It should be considered to start with the smallest size Distractor to avoid over distraction. Apply the Distractor parallel to the intervertebral space and turn clockwise to open up the disc space.



Figure 4
Correct the surface

To obtain parallel distraction, proceed in the counter side with one size bigger until sufficient distraction is achieved. Correct the surface if necessary with the Rasp tool (Figure 4).



Figure 5
Determination the implant size

Determination of Implant Size

The selection of the implant trial depends on the height, width and depth of the intervertebral disc space, as well as the chosen preparation technique and the anatomy of the patient (Figure 5). Select the trial that best matches the shape of the intervertebral disc space. Choose gradually the trial that is higher and/or with higher lordosis angle until a secure position in the intervertebral disc space is achieved. Use the highest possible trial implant to ensure maximal stability. You can verify the correct position of the implant trial with help of the x-ray imaging to make sure the implant fits exactly into the intervertebral disc space. The length of the implant can be determined according to the picture below (Figure 5a).




Figure 5a



Figure 6
Implant mounting

Implant Insertion

Open the sterile packaging of the implant size that was determined with the Trial. Attach the implant to the Inserter. Place the implant between the tips of the Inserter that fit into the notches on the lateral side of the implant. Screw the implant to the Inserter by final tightening of the screwcap on the back of the Inserter and tighten with hand (Figure 6). Once the implant is blocked on the instrument, it can be implanted into the intervertebral space.



Figure 7
Place the implant

It is recommended to place the implants in pairs (Figure 7). In case of implanting an oblique single implant, it is advisable to use the largest possible footprint available that fits the anatomy.

Gently tap the implant into the intervertebral disc space using the hammer provided in the instrument set. The first implant should not cross the midline in order not to interfere with the second implant.

During implantation it is necessary to retract and protect neural structures by a nerve root retractor. Despite the application of smooth solid titanium surfaces, the rough porous structure potentially increases the risk of soft tissue adherence unless there is sufficient protection. Use the Hexagonal Key to rotate the implant if it needs to be assembler (Figure 8).



Figure 8
Implant rotation

Completion of Surgery

After Implantation of the implant a final check of the position of the implant under fluososcopy is advised (Figure 9).

Posterior supplemental fixation is strongly recommended. Perform wound closure as usual. Please document which implants were used in the patient files. Patient labels are supplied with each implant for your convenience. The use of the BMD Titanium PLIF GW Implant[®] does not require any specific postoperative care and the patient should be treated according to hospital and medical standards.




Figure 9
Final position



Figure 10
Implant removal

Implant Removal

Either the Inserter or Universal Removal Instrument may be used for Implant removal by attachment via clockwise rotation to the implant threads (Figure 10). Be careful to avoid pushing the implant posteriorly. Once the implant is firmly attached, remove the implant from the disc space. Vertebral bone overgrowth or osteophytes may be removed to facilitate implant retrieval.

List of Implants

DESCRIPTION	L x W x H	ANGLE	CODE
PLIF GW	33mm x 12mm x 7mm	0°	BMDH33120700
PLIF GW	33mm x 12mm x 8mm	0°	BMDH33120800
PLIF GW	33mm x 12mm x 9mm	0°	BMDH33120900
PLIF GW	33mm x 12mm x 10mm	0°	BMDH33121000
PLIF GW	33mm x 12mm x 11mm	0°	BMDH33121100
PLIF GW	33mm x 12mm x 12mm	0°	BMDH33121200
PLIF GW	33mm x 12mm x 13mm	0°	BMDH33121300
PLIF GW	33mm x 12mm x 14mm	0°	BMDH33121400
PLIF GW	33mm x 12mm x 15mm	0°	BMDH33121500
PLIF GW	20mm x 9mm x 7mm	4°	BMDH20090704
PLIF GW	20mm x 9mm x 8mm	4°	BMDH20090804
PLIF GW	20mm x 9mm x 9mm	4°	BMDH20090904
PLIF GW	20mm x 9mm x 10mm	4°	BMDH20091004
PLIF GW	20mm x 9mm x 11mm	4°	BMDH20091104
PLIF GW	20mm x 9mm x 12mm	4°	BMDH20091204
PLIF GW	20mm x 9mm x 13mm	4°	BMDH20091304
PLIF GW	20mm x 9mm x 14mm	4°	BMDH20091404
PLIF GW	20mm x 9mm x 15mm	4°	BMDH20091504
PLIF GW	20mm x 9mm x 16mm	4°	BMDH20091604
PLIF GW	20mm x 9mm x 17mm	4°	BMDH20091704
PLIF GW	24mm x 9mm x 7mm	4°	BMDH24090704
PLIF GW	24mm x 9mm x 8mm	4°	BMDH24090804
PLIF GW	24mm x 9mm x 9mm	4°	BMDH24090904
PLIF GW	24mm x 9mm x 10mm	4°	BMDH24091004
PLIF GW	24mm x 9mm x 11mm	4°	BMDH24091104
PLIF GW	24mm x 9mm x 12mm	4°	BMDH24091204
PLIF GW	24mm x 9mm x 13mm	4°	BMDH24091304
PLIF GW	24mm x 9mm x 14mm	4°	BMDH24091404
PLIF GW	24mm x 9mm x 15mm	4°	BMDH24091504
PLIF GW	24mm x 9mm x 16mm	4°	BMDH24091604
PLIF GW	24mm x 9mm x 17mm	4°	BMDH24091704
PLIF GW	24mm x 10mm x 7mm	4°	BMDH24100704
PLIF GW	24mm x 10mm x 8mm	4°	BMDH24100804
PLIF GW	24mm x 10mm x 9mm	4°	BMDH24100904
PLIF GW	24mm x 10mm x 10mm	4°	BMDH24101004
PLIF GW	24mm x 10mm x 11mm	4°	BMDH24101104
PLIF GW	24mm x 10mm x 12mm	4°	BMDH24101204
PLIF GW	24mm x 10mm x 13mm	4°	BMDH24101304
PLIF GW	24mm x 10mm x 14mm	4°	BMDH24101404
PLIF GW	24mm x 10mm x 15mm	4°	BMDH24101504
PLIF GW	24mm x 10mm x 16mm	4°	BMDH24101604
PLIF GW	24mm x 10mm x 17mm	4°	BMDH24101704
PLIF GW	28mm x 9mm x 7mm	4°	BMDH28090704
PLIF GW	28mm x 9mm x 8mm	4°	BMDH28090804
PLIF GW	28mm x 9mm x 9mm	4°	BMDH28090904
PLIF GW	28mm x 9mm x 10mm	4°	BMDH28091004
PLIF GW	28mm x 9mm x 11mm	4°	BMDH28091104
PLIF GW	28mm x 9mm x 12mm	4°	BMDH28091204

DESCRIPTION	L x W x H	ANGLE	CODE
PLIF GW	28mm x 9mm x 13mm	4°	BMDH28091304
PLIF GW	28mm x 9mm x 14mm	4°	BMDH28091404
PLIF GW	28mm x 9mm x 15mm	4°	BMDH28091504
PLIF GW	28mm x 9mm x 16mm	4°	BMDH28091604
PLIF GW	28mm x 9mm x 17mm	4°	BMDH28091704
PLIF GW	28mm x 10mm x 7mm	4°	BMDH28100704
PLIF GW	28mm x 10mm x 8mm	4°	BMDH28100804
PLIF GW	28mm x 10mm x 9mm	4°	BMDH28100904
PLIF GW	28mm x 10mm x 10mm	4°	BMDH28101004
PLIF GW	28mm x 10mm x 11mm	4°	BMDH28101104
PLIF GW	28mm x 10mm x 12mm	4°	BMDH28101204
PLIF GW	28mm x 10mm x 13mm	4°	BMDH28101304
PLIF GW	28mm x 10mm x 14mm	4°	BMDH28101404
PLIF GW	28mm x 10mm x 15mm	4°	BMDH28101504
PLIF GW	28mm x 10mm x 16mm	4°	BMDH28101604
PLIF GW	28mm x 10mm x 17mm	4°	BMDH28101704
PLIF GW	33mm x 12mm x 9mm	4°	BMDH33120904
PLIF GW	33mm x 12mm x 11mm	4°	BMDH33121104
PLIF GW	33mm x 12mm x 13mm	4°	BMDH33121304
PLIF GW	33mm x 12mm x 15mm	4°	BMDH33121504
PLIF GW	20mm x 9mm x 7mm	8°	BMDH20090708
PLIF GW	20mm x 9mm x 8mm	8°	BMDH20090808
PLIF GW	20mm x 9mm x 9mm	8°	BMDH20090908
PLIF GW	20mm x 9mm x 10mm	8°	BMDH20091008
PLIF GW	20mm x 9mm x 11mm	8°	BMDH20091108
PLIF GW	20mm x 9mm x 12mm	8°	BMDH20091208
PLIF GW	20mm x 9mm x 13mm	8°	BMDH20091308
PLIF GW	20mm x 9mm x 14mm	8°	BMDH20091408
PLIF GW	20mm x 9mm x 15mm	8°	BMDH20091508
PLIF GW	20mm x 9mm x 16mm	8°	BMDH20091608
PLIF GW	20mm x 9mm x 17mm	8°	BMDH20091708
PLIF GW	24mm x 9mm x 7mm	8°	BMDH24090708
PLIF GW	24mm x 9mm x 8mm	8°	BMDH24090808
PLIF GW	24mm x 9mm x 9mm	8°	BMDH24090908
PLIF GW	24mm x 9mm x 10mm	8°	BMDH24091008
PLIF GW	24mm x 9mm x 11mm	8°	BMDH24091108
PLIF GW	24mm x 9mm x 12mm	8°	BMDH24091208
PLIF GW	24mm x 9mm x 13mm	8°	BMDH24091308
PLIF GW	24mm x 9mm x 14mm	8°	BMDH24091408
PLIF GW	24mm x 9mm x 15mm	8°	BMDH24091508
PLIF GW	24mm x 9mm x 16mm	8°	BMDH24091608
PLIF GW	24mm x 9mm x 17mm	8°	BMDH24091708
PLIF GW	24mm x 10mm x 7mm	8°	BMDH24100708
PLIF GW	24mm x 10mm x 8mm	8°	BMDH24100808
PLIF GW	24mm x 10mm x 9mm	8°	BMDH24100908
PLIF GW	24mm x 10mm x 10mm	8°	BMDH24101008
PLIF GW	24mm x 10mm x 11mm	8°	BMDH24101108
PLIF GW	24mm x 10mm x 12mm	8°	BMDH24101208

DESCRIPTION	L x W x H	ANGLE	CODE
PLIF GW	24mm x 10mm x 13mm	8°	BMDH24101308
PLIF GW	24mm x 10mm x 14mm	8°	BMDH24101408
PLIF GW	24mm x 10mm x 15mm	8°	BMDH24101508
PLIF GW	24mm x 10mm x 16mm	8°	BMDH24101608
PLIF GW	24mm x 10mm x 17mm	8°	BMDH24101708
PLIF GW	28mm x 9mm x 7mm	8°	BMDH28090708
PLIF GW	28mm x 9mm x 8mm	8°	BMDH28090808
PLIF GW	28mm x 9mm x 9mm	8°	BMDH28090908
PLIF GW	28mm x 9mm x 10mm	8°	BMDH28091008
PLIF GW	28mm x 9mm x 11mm	8°	BMDH28091108
PLIF GW	28mm x 9mm x 12mm	8°	BMDH28091208
PLIF GW	28mm x 9mm x 13mm	8°	BMDH28091308
PLIF GW	28mm x 9mm x 14mm	8°	BMDH28091408
PLIF GW	28mm x 9mm x 15mm	8°	BMDH28091508
PLIF GW	28mm x 9mm x 16mm	8°	BMDH28091608
PLIF GW	28mm x 9mm x 17mm	8°	BMDH28091708
PLIF GW	28mm x 10mm x 7mm	8°	BMDH28100708
PLIF GW	28mm x 10mm x 8mm	8°	BMDH28100808
PLIF GW	28mm x 10mm x 9mm	8°	BMDH28100908
PLIF GW	28mm x 10mm x 10mm	8°	BMDH28101008
PLIF GW	28mm x 10mm x 11mm	8°	BMDH28101108
PLIF GW	28mm x 10mm x 12mm	8°	BMDH28101208
PLIF GW	28mm x 10mm x 13mm	8°	BMDH28101308
PLIF GW	28mm x 10mm x 14mm	8°	BMDH28101408
PLIF GW	28mm x 10mm x 15mm	8°	BMDH28101508
PLIF GW	28mm x 10mm x 16mm	8°	BMDH28101608
PLIF GW	28mm x 10mm x 17mm	8°	BMDH28101708
PLIF GW	33mm x 12mm x 9mm	8°	BMDH33120908
PLIF GW	33mm x 12mm x 10mm	8°	BMDH33121008
PLIF GW	33mm x 12mm x 11mm	8°	BMDH33121108
PLIF GW	33mm x 12mm x 12mm	8°	BMDH33121208
PLIF GW	33mm x 12mm x 13mm	8°	BMDH33121308
PLIF GW	33mm x 12mm x 14mm	8°	BMDH33121408
PLIF GW	33mm x 12mm x 15mm	8°	BMDH33121508
PLIF GW	20mm x 9mm x 7mm	12°	BMDH20090712
PLIF GW	20mm x 9mm x 8mm	12°	BMDH20090812
PLIF GW	20mm x 9mm x 9mm	12°	BMDH20090912
PLIF GW	20mm x 9mm x 10mm	12°	BMDH20091012
PLIF GW	20mm x 9mm x 11mm	12°	BMDH20091112
PLIF GW	20mm x 9mm x 12mm	12°	BMDH20091212
PLIF GW	20mm x 9mm x 13mm	12°	BMDH20091312
PLIF GW	20mm x 9mm x 14mm	12°	BMDH20091412
PLIF GW	20mm x 9mm x 15mm	12°	BMDH20091512
PLIF GW	20mm x 9mm x 16mm	12°	BMDH20091612
PLIF GW	20mm x 9mm x 17mm	12°	BMDH20091712
PLIF GW	24mm x 9mm x 9mm	12°	BMDH24090912
PLIF GW	24mm x 9mm x 10mm	12°	BMDH24091012
PLIF GW	24mm x 9mm x 11mm	12°	BMDH24091112
PLIF GW	24mm x 9mm x 12mm	12°	BMDH24091212
PLIF GW	24mm x 9mm x 13mm	12°	BMDH24091312

DESCRIPTION	L x W x H	ANGLE	CODE
PLIF GW	24mm x 9mm x 14mm	12°	BMDH24091412
PLIF GW	24mm x 9mm x 15mm	12°	BMDH24091512
PLIF GW	24mm x 9mm x 16mm	12°	BMDH24091612
PLIF GW	24mm x 9mm x 17mm	12°	BMDH24091712
PLIF GW	24mm x 10mm x 9mm	12°	BMDH24100912
PLIF GW	24mm x 10mm x 10mm	12°	BMDH24101012
PLIF GW	24mm x 10mm x 11mm	12°	BMDH24101112
PLIF GW	24mm x 10mm x 12mm	12°	BMDH24101212
PLIF GW	24mm x 10mm x 13mm	12°	BMDH24101312
PLIF GW	24mm x 10mm x 14mm	12°	BMDH24101412
PLIF GW	24mm x 10mm x 15mm	12°	BMDH24101512
PLIF GW	24mm x 10mm x 16mm	12°	BMDH24101612
PLIF GW	24mm x 10mm x 17mm	12°	BMDH24101712
PLIF GW	28mm x 9mm x 11mm	12°	BMDH28091112
PLIF GW	28mm x 9mm x 12mm	12°	BMDH28091212
PLIF GW	28mm x 9mm x 13mm	12°	BMDH28091312
PLIF GW	28mm x 9mm x 14mm	12°	BMDH28091412
PLIF GW	28mm x 9mm x 15mm	12°	BMDH28091512
PLIF GW	28mm x 9mm x 16mm	12°	BMDH28091612
PLIF GW	28mm x 9mm x 17mm	12°	BMDH28091712
PLIF GW	28mm x 10mm x 11mm	12°	BMDH28101112
PLIF GW	28mm x 10mm x 12mm	12°	BMDH28101212
PLIF GW	28mm x 10mm x 13mm	12°	BMDH28101312
PLIF GW	28mm x 10mm x 14mm	12°	BMDH28101412
PLIF GW	28mm x 10mm x 15mm	12°	BMDH28101512
PLIF GW	28mm x 10mm x 16mm	12°	BMDH28101612
PLIF GW	28mm x 10mm x 17mm	12°	BMDH28101712
PLIF GW	33mm x 12mm x 11mm	12°	BMDH33121112
PLIF GW	33mm x 12mm x 12mm	12°	BMDH33121212
PLIF GW	33mm x 12mm x 13mm	12°	BMDH33121312
PLIF GW	33mm x 12mm x 14mm	12°	BMDH33121412
PLIF GW	33mm x 12mm x 15mm	12°	BMDH33121512
PLIF GW	33mm x 12mm x 16mm	12°	BMDH33121612
PLIF GW	33mm x 12mm x 17mm	12°	BMDH33121712
PLIF GW	20mm x 9mm x 7mm	16°	BMDH20090716
PLIF GW	20mm x 9mm x 8mm	16°	BMDH20090816
PLIF GW	20mm x 9mm x 9mm	16°	BMDH20090916
PLIF GW	20mm x 9mm x 10mm	16°	BMDH20091016
PLIF GW	20mm x 9mm x 11mm	16°	BMDH20091116
PLIF GW	20mm x 9mm x 12mm	16°	BMDH20091216
PLIF GW	20mm x 9mm x 13mm	16°	BMDH20091316
PLIF GW	20mm x 9mm x 14mm	16°	BMDH20091416
PLIF GW	20mm x 9mm x 15mm	16°	BMDH20091516
PLIF GW	20mm x 9mm x 16mm	16°	BMDH20091616
PLIF GW	20mm x 9mm x 17mm	16°	BMDH20091716
PLIF GW	24mm x 9mm x 11mm	16°	BMDH24091116
PLIF GW	24mm x 9mm x 12mm	16°	BMDH24091216
PLIF GW	24mm x 9mm x 13mm	16°	BMDH24091316
PLIF GW	24mm x 9mm x 14mm	16°	BMDH24091416
PLIF GW	24mm x 9mm x 15mm	16°	BMDH24091516

DESCRIPTION	L x W x H	ANGLE	CODE
PLIF GW	24mm x 9mm x 16mm	16°	BMDH24091616
PLIF GW	24mm x 9mm x 17mm	16°	BMDH24091716
PLIF GW	24mm x 10mm x 11mm	16°	BMDH24101116
PLIF GW	24mm x 10mm x 12mm	16°	BMDH24101216
PLIF GW	24mm x 10mm x 13mm	16°	BMDH24101316
PLIF GW	24mm x 10mm x 14mm	16°	BMDH24101416
PLIF GW	24mm x 10mm x 15mm	16°	BMDH24101516
PLIF GW	24mm x 10mm x 16mm	16°	BMDH24101616
PLIF GW	24mm x 10mm x 17mm	16°	BMDH24101716
PLIF GW	28mm x 9mm x 13mm	16°	BMDH28091316
PLIF GW	28mm x 9mm x 14mm	16°	BMDH28091416
PLIF GW	28mm x 9mm x 15mm	16°	BMDH28091516
PLIF GW	28mm x 9mm x 16mm	16°	BMDH28091616
PLIF GW	28mm x 9mm x 17mm	16°	BMDH28091716
PLIF GW	28mm x 10mm x 13mm	16°	BMDH28101316
PLIF GW	28mm x 10mm x 14mm	16°	BMDH28101416
PLIF GW	28mm x 10mm x 15mm	16°	BMDH28101516
PLIF GW	28mm x 10mm x 16mm	16°	BMDH28101616
PLIF GW	28mm x 10mm x 17mm	16°	BMDH28101716
PLIF GW	20mm x 9mm x 9mm	20°	BMDH20090920
PLIF GW	20mm x 9mm x 10mm	20°	BMDH20091020
PLIF GW	20mm x 9mm x 11mm	20°	BMDH20091120
PLIF GW	20mm x 9mm x 12mm	20°	BMDH20091220
PLIF GW	20mm x 9mm x 13mm	20°	BMDH20091320
PLIF GW	20mm x 9mm x 14mm	20°	BMDH20091420
PLIF GW	20mm x 9mm x 15mm	20°	BMDH20091520
PLIF GW	20mm x 9mm x 16mm	20°	BMDH20091620
PLIF GW	20mm x 9mm x 17mm	20°	BMDH20091720
PLIF GW	24mm x 9mm x 13mm	20°	BMDH24091320
PLIF GW	24mm x 9mm x 14mm	20°	BMDH24091420
PLIF GW	24mm x 9mm x 15mm	20°	BMDH24091520
PLIF GW	24mm x 9mm x 16mm	20°	BMDH24091620
PLIF GW	24mm x 9mm x 17mm	20°	BMDH24091720
PLIF GW	24mm x 10mm x 13mm	20°	BMDH24101320
PLIF GW	24mm x 10mm x 14mm	20°	BMDH24101420
PLIF GW	24mm x 10mm x 15mm	20°	BMDH24101520
PLIF GW	24mm x 10mm x 16mm	20°	BMDH24101620
PLIF GW	24mm x 10mm x 17mm	20°	BMDH24101720
PLIF GW	28mm x 9mm x 15mm	20°	BMDH28091520
PLIF GW	28mm x 9mm x 16mm	20°	BMDH28091620
PLIF GW	28mm x 9mm x 17mm	20°	BMDH28091720
PLIF GW	28mm x 10mm x 15mm	20°	BMDH28101520
PLIF GW	28mm x 10mm x 16mm	20°	BMDH28101620
PLIF GW	28mm x 10mm x 17mm	20°	BMDH28101720

Tool List

NAME	REF	NAME	REF
Curette D9L5	35-110	PLIF Trial 0°, 13 mm	45-756
Curette 45° D9L5	35-111	PLIF Trial 0°, 14 mm	45-757
Curette 15° D7L3	35-112	PLIF Trial 0°, 15 mm	45-758
Curette 15° D9L5	35-113	PLIF Trial 4°, 7 mm	45-759
Cage Holder	45-103	PLIF Trial 4°, 8 mm	45-760
Mallet	35-115	PLIF Trial 4°, 9 mm	45-761
Rasp	35-116	PLIF Trial 4°, 10 mm	45-762
Removal Instrument	35-119	PLIF Trial 4°, 11 mm	45-763
T-Handle	35-114	PLIF Trial 4°, 12 mm	45-764
Final Impactor	45-101	PLIF Trial 4°, 13 mm	45-765
Hexagonal Key	45-102	PLIF Trial 4°, 14 mm	45-766
Cage Slider	45-104	PLIF Trial 4°, 15 mm	45-767
Distractor 7 mm	35-130	PLIF Trial 4°, 16 mm	45-768
Distractor 8 mm	35-131	PLIF Trial 4°, 17 mm	45-769
Distractor 9 mm	35-132	PLIF Trial 8°, 7 mm	45-770
Distractor 10 mm	35-133	PLIF Trial 8°, 8 mm	45-771
Distractor 11 mm	35-134	PLIF Trial 8°, 9 mm	45-772
Distractor 12 mm	35-135	PLIF Trial 8°, 10 mm	45-773
Distractor 13 mm	35-136	PLIF Trial 8°, 11 mm	45-774
Distractor 14 mm	35-137	PLIF Trial 8°, 12 mm	45-775
Distractor 15 mm	35-138	PLIF Trial 8°, 13 mm	45-776
Distractor 16 mm	35-139	PLIF Trial 8°, 14 mm	45-777
Distractor 17 mm	35-140	PLIF Trial 8°, 15 mm	45-778
Distractor shaver 7 mm	35-101	PLIF Trial 8°, 16 mm	45-779
Distractor shaver 8 mm	35-102	PLIF Trial 8°, 17 mm	45-780
Distractor shaver 9 mm	35-103	PLIF Trial 12°, 9 mm	45-781
Distractor shaver 10 mm	35-104	PLIF Trial 12°, 10 mm	45-782
Distractor shaver 11 mm	35-105	PLIF Trial 12°, 11 mm	45-783
Distractor shaver 12 mm	35-106	PLIF Trial 12°, 12 mm	45-784
Distractor shaver 13 mm	35-107	PLIF Trial 12°, 13 mm	45-785
Distractor shaver 14 mm	35-108	PLIF Trial 12°, 14 mm	45-786
Distractor shaver 15 mm	35-109	PLIF Trial 12°, 15 mm	45-787
Distractor shaver 16 mm	35-117	PLIF Trial 12°, 16 mm	45-788
Distractor shaver 17 mm	35-118	PLIF Trial 12°, 17 mm	45-789
Shaver 7 mm	35-150	PLIF Trial 16°, 11 mm	45-790
Shaver 8 mm	35-151	PLIF Trial 16°, 12 mm	45-791
Shaver 9 mm	35-152	PLIF Trial 16°, 13 mm	45-792
Shaver 10 mm	35-153	PLIF Trial 16°, 14 mm	45-793
Shaver 11 mm	35-154	PLIF Trial 16°, 15 mm	45-794
Shaver 12 mm	35-155	PLIF Trial 16°, 16 mm	45-795
Shaver 13 mm	35-156	PLIF Trial 16°, 17 mm	45-796
Shaver 14 mm	35-157	PLIF Trial 20°, 13 mm	45-797
Shaver 15 mm	35-158	PLIF Trial 20°, 14 mm	45-798
Shaver 16 mm	35-159	PLIF Trial 20°, 15 mm	45-799
Shaver 17 mm	35-160	PLIF Trial 20°, 16 mm	45-800
PLIF Trial 0°, 7 mm	45-750	PLIF Trial 20°, 17 mm	45-801
PLIF Trial 0°, 8 mm	45-751	Trial Basket	35-120
PLIF Trial 0°, 9 mm	45-752	Instrument Basket	45-400
PLIF Trial 0°, 10 mm	45-753		
PLIF Trial 0°, 11 mm	45-754		
PLIF Trial 0°, 12 mm	45-755		

Important Information on the Global BMD Titanium PLIF GW Implant[®]

Device Description and Materials

Global Biomedica offers spinal surgery implants with excellent biocompatibility and bioactivity. We use the latest innovative technologies to create a line of spinal BMD Titanium Implant[®]. These are implant for interbody fusion of biocompatible titanium (Ti64ELI). Internal and surface grid structure with optimum pore size of 700µm (Cube vertex centroid - lattice) with reinforced edges ensures not only high stability and resistance to deformation of the implant, or immersion in the vertebral body, but also bioactivity - potentiates the formation of bone in the area of contact surfaces and thus the formation of a strong connection between the implant and bone, the risk of developing non union (pseudoarthrosis) is thus minimized.

Indication for Use

The lumbar BMD Titanium Implant[®] is indicated for lumbar interbody fusion for the following indications:

- Degenerative and postoperative spinal involvement in term of rupture, osteochondrosis, degenerative stenosis, pseudarthrosis
- Traumatic disc involvement – only in the absence of traumatic vertebrae involvement and significant ligament instability of posterior column

Contraindications

The lumbar BMD Titanium Implant[®] should not be used in patients with any of the following contraindications (Contraindications include, but are not limited). Infection or progressive infection, fever or inflammation.

- Obesity
- Insanity
- Allergy for any of system components
- Any anatomical, medical or surgical conditions that may prevent the potential or intentional benefits of using spinal implants
- Bone, joint or ligament conditions such as: Osteolytic involvement, bone absorption, osteomalacia. Osteopenia and osteoporosis is a relative contraindication and a thorough assessment of the situation is required before performing surgery
- Implant size, shape or functionality of the anchor may not be sufficient to achieve the expected clinical results
- Any risk of patient unwillingness to follow postoperative instructions
- Any other than described indications

Warnings And Precautions

- The effectiveness and safety of interbody fusion is only applicable for certain conditions with significant instability which require the fusion supported by medical device.
- Correct placement and appropriate size selection are crucial to achieve optimal results.
- The device might be supportive for such mechanical instability like deformity, fracture, listhesis, dislocation, tumor, pseudoarthrosis. The effectiveness and safety for any other conditions are unknown.
- BMD Titanium Implant[®] may be supported by additional fixation device. In some cases additional fixation device is highly recommended.
- The contemporary applications of pedicular fixation must be performed by experienced surgeons with specific training in usage BMD Titanium Implant[®].
- The spinal screw fixation system and/or interbody cage system should not be considered as sole spinal support. No implants can withstand body loads without bone support.
- Over the time bends, breakages, loosening, migration may occur. Successful results are not always achievable. The factors as proper pre-operative and operative procedure, comprehensive knowledge or surgical techniques, proper selection of implant's size and type are considerably important in treatment process.
- Patients with obesity, smokers, alcohol abused are enhanced risk of non-fusion. Also patients in weak muscle or bone conditions, nervous system dysfunctions are poor candidates for spinal fusion.
- Prior or during or after the surgery in order to evaluate or check positioning of the implants or patients anatomy or any other patients or implants X-ray or CT or any other invasive diagnostic examinations may be necessary to be performed.
- The proper, patient's individual implants selection in terms of type, size, shape or instruments handling is crucial.
- Extensive bending or contouring should be avoided.
- Sharp edges cutting, reversed bending, scratching or notching may generate internal stressing which may weak the implants or construct.

Potential Adverse Effects

Potential adverse events which may occur after the spinal surgery with or without instrumentation include, but are not limited:

- Breaking, bending, and/or breakage of any or all of the system components.
- Migrations of any system components.
- Pressure on the skin from component parts in patients with inadequate tissue coverage.
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments.
- Damage and leakage of the hard arachnoid, leakage of coeliolymph.
- Neurological dysfunctions such as paresthesia, radiculopathy, paralysis, hyperesthesia or any others related to general surgery associated to anesthesia.
- Infections.
- Loss of urinary and defaecate functions.
- Permanent or temporary or developing sexual dysfunctions.
- Postoperative change in body curvature, change of physiological range of movement.
- Pseudarthrosis or non-fusion or delayed fusion.
- Loss of bone or overgrowth.
- Permanent or temporary limitation or inability to perform daily activities.
- Changing in mental behavior.
- Permanent or temporary or developing respiratory problems.
- Permanent or temporary or developing cardiovascular deteriorations or dysfunctions.
- Death.

In some cases additional surgery or surgeries might be necessary to correct or change potential adverse events.

IMPORTANT NOTICE: All necessary information on surgery, potential risks, benefits and adverse effects must be reported to the patient prior to surgery.



GLOBAL[®]
biomedica

Disclaimer: This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is recommendations. Because this information does not purport to constitute any diagnostics or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part.



Caution: Please see the product Instructions for Use for a complete listing of the indications, contraindications, precautions, warnings and adverse effects.



MANUFACTURER
GLOBAL BIOMEDICA s.r.o., Ostravská 555/24,
737 01 Český Těšín, Czech Republic, info@globalbmd.com

©2020 GLOBAL BIOMEDICA s.r.o. All rights reserved.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to GLOBAL BIOMEDICA or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of GLOBAL BIOMEDICA. This material is intended for health care professionals, the GLOBAL BIOMEDICA sales force and authorized representatives. Distribution to any other recipient is prohibited. For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert and www.globalbmd.com. Check for country product clearances and reference product specific instructions for use.

Surgical Technique Guide – BMD Titanium PLIF GW Implant[®], Rev. 2021-04-21 v. 2.1