
Erisma[®]-LP

Surgical technique

Inspired by the spine

Inspired for the spine

Contents

Introduction	3
Features & benefits	4
Surgical technique	6
Product catalog	14
Instructions for use	18

Introduction

Designed by a team with extensive expertise in the development of spinal instrumentation, Erisma–LP is the result of optimized development based on an analysis of existing technologies.

The Erisma–LP spinal instrumentation is a pedicle screw based fixation. The design team focused on creating a top loading system that incorporates reduced implant volume and ergonomically designed instrumentation. Erisma–LP was developed to meet the spine surgeon’s needs for reliability, security and ease of use during day-to-day surgical practice. Erisma–LP is a

comprehensive system of implants and instruments designed to address pathologies of the spine in the thoraco-lumbar region. Composed of polyaxial screws, monoaxial screws, crosslinks, rods and set screws, Erisma–LP integrates a wide selection of sizes and diameters, allowing the surgeon to perform posterior fixation of the spine for various pathologies.

› Indications:

Disc Degenerative disease
 Spondylolisthesis
 Spinal stenosis
 Lumbar deformities
 Trauma
 Tumor
 Pseudarthrosis
 Failed previous fusion


Features & benefits

SCREWS

Optimized volume and profile

Limited interference with anatomical features

Tapered thread

Enhanced strength under bending moments

Self-tapping blunt tip

Easy, smooth and non-aggressive screw insertion

Color-coded screws and taps


Quick screw selection

Buttress thread

No splaying of the screw head

Optimized performance of the screw/rod connection

Optimized thread entry

Prevents cross-threading


SET SCREW

RODS

Short pre-bent rods

Quick rod preparation

Straight rods with guiding line

Facilitate bending of the rod within one plane

Facilitate rod positioning

CROSSLINK

Fully adjustable crosslink

Easy implantation in all configurations

INSTRUMENTATION

User-friendly

Precise, quick and efficient maneuvers


Quick connect and modular handles

Adaptability to the surgeon's habits

Ergonomic design

Complete instrumentation, compact instrument set

Surgical Technique

- | | | | |
|---|--------------------------------------|----|-----------------------|
| 1 | Pedicle preparation | 7 | Set screw positioning |
| 2 | Pedicle drilling and screw selection | 8 | Correction maneuvers |
| 3 | Tapping | 9 | Final tightening |
| 4 | Pedicle screw insertion | 10 | Crosslink placement |
| 5 | Rod preparation | 11 | Removal |
| 6 | Rod insertion | | |

Following all necessary safety protocols, the patient is positioned on the operating table in the prone position. The surgical approach is performed according to the known and validated techniques routinely used by surgeons.


Pedicle preparation

1

99702000 SQUARE AWL


After adequate exposure, the appropriate pedicle entry point is located and the entrance to the pedicle canal is perforated using the **square awl**.

Pedicle drilling and screw selection


2

A pathway in the intrapedicular cancellous bone is then opened up using the **pedicle probe** at the appropriate angle and depth. The pedicle probe is graduated with 10mm intervals to determine the depth measurement, as well as to determine the proper screw length.

99731000 PEDICLE PROBE

The integrity of the pedicular hole's walls is then checked with the **pedicle sounder**.



99730000 PEDICLE SOUNDER

Tapping

3

The Erisma-LP screws are self-tapping and do not require tapping. For the surgeon's convenience, optional tapping can be performed using the **taps** provided in all diameters. The screw diameter is determined by the color-code indicated on the tap.

997010** TAP

18710002
HOLDING SCREWDRIVER

HANDLES

Pedicle screw insertion

4



The appropriate screw diameter and length are determined.

i The proper screw length may be checked using the screw gauge placed on the instrument tray.

To implant the screw, fully engage the tip of the **holding screwdriver** in the head of the screw ① and turn the locking knob clockwise ②. Then, make sure that the screw is straight and secured ③.


The screw is inserted in the pedicle hole and the holding screwdriver is disengaged by turning the locking knob counterclockwise. Repeat for each screw.

Rod preparation

Once the screws have been inserted and their positions have been verified by X-ray, the appropriate rod length is selected using the **rod template**. If necessary, bend the rod using the **French bender** to set the appropriate spinal contours.


Rod insertion

The rod is inserted in the screw heads using the **angled rod holding forceps**. The Erisma-LP system offers two options for linking the rod to the screw heads:

› Rod Pusher

When the rod is slightly above the screw head, the **rod pusher** may be used to push the rod into position.

99733*** ROD TEMPLATE

99750001 FRENCH BENDER

99715003
NARROW NOSE ROD HOLDER

99715002
ANGLED ROD HOLDING FORCEPS


99743000 ROD PUSHER

18742002 PERSUADER

► Persuader


The **persuader** is used if additional force is needed to bring the rod to the implant. The persuader is connected to the head of the screw. The trigger is then pressed and the tube slid down.

Once the lever is pushed down, the rod is fully seated in the implant. When the introduction lever is tipped down, the persuader stays in position.

i *The lever can be pushed down only when the persuader is well-positioned on the implant and the tube slid down.*

When finished, press the trigger again, slide the tube up and pull the lever up. The persuader can be removed and is ready to use again.


01714001 SET SCREW HOLDER

Set screw positioning

The set screws are fitted on the tip of the **set screw holder** and sequentially threaded on the implant. Only a few threads are necessary to secure the rod in the implants. Repeat for each set screw.

i *It is recommended to use the persuader to guarantee that the rod is perfectly in position and therefore, avoid possible cross-threading.*


Correction maneuvers

8


Before final tightening of the construct, rod rotation, compression or distraction may be achieved.

› Rod rotation

To achieve the final rod positioning, the rod can be turned using the [narrow nose rod holder](#).


99715003
NARROW NOSE ROD HOLDER

99741001 DISTRACTOR

› Distraction and compression

Segmental distraction or compression may be carried out using the [compressor](#) or the [distractor](#).

99740001 COMPRESSOR

01721003 COUNTER TORQUE


99781200
TORQUE LIMITED T HANDLE01720003
TIGHTENING WRENCH T30

18733001 CROSSLINK TEMPLATE

Final tightening

9


The set screw is tightened using the **tightening wrench T30** attached to the **torque limited T-handle**. The **counter torque** is positioned over the screw head and the wrench is inserted through the tube of the counter-torque. Tighten until the torque limited T-handle triggers once. Repeat for each screw.


Crosslink placement

10




To complete the construct, a crosslink may be added to provide increased stability. The appropriate size of the crosslink (S, M, L, XL, XXL) is determined using the **crosslink template**.

Once the crosslink size is determined and its unlocked position is checked, the **crosslink holder** is used to position the crosslink. The first transverse connector is positioned and its set screw is tightened using the **tightening wrench** ①. The operation is repeated for the second connector ②. When lateral set screws are secured, the center set screw is tightened ③.

18713001 CROSSLINK HOLDER

01720001
TIGHTENING WRENCH T20

Removal

11

The set screws can be removed using the **counter torque** and the **tightening wrench** to turn the set screws counterclockwise. When the set screws are removed, the rod can be easily extracted using the **angled rod holding forceps**. The **holding screwdriver** is used to remove the screws.

Product catalog

SQUARE AWL

99702000

PEDICLE SOUNDER

99730000

PEDICLE PROBE

99731000

TAP

Ø4.5mm	99701045
Ø5.5mm	99701055
Ø6.5mm	99701065
Ø7.5mm	99701075
Ø8.5mm	99701085

HOLDING SCREWDRIVER

18710002

**FIXED PALM HANDLE
NON RATCHETING T-HANDLE
RATCHETING HANDLE
TORQUE LIMITED T-HANDLE**

99780000
99781000
99781100
99781200

ANGLED ROD HOLDING FORCEPS

99715002

FRENCH BENDER

99750001

NARROW NOSE ROD HOLDER

99715003

ROD PUSHER

99743000

PERSUADER

18742002

SET SCREW HOLDER

01714001

DISTRACTOR

99741001

COMPRESSOR

99740001


COUNTER TORQUE

01721003

TIGHTENING WRENCH T30
TIGHTENING WRENCH T20

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



EN Instructions for use

The Erisma-LP instrumentation is designed for the surgical treatment of spinal pathologies. The treatment consists of the fusion of two or several vertebrae in order to restore spinal stability, with or without any other endocanalicular concomitant surgical procedure. For optimal results, a detailed preoperative evaluation, a meticulous surgical technique and adequate post-operative care are mandatory. It is important that both the patient and surgeon are fully aware of the risks and possible complications associated with this type of surgery. This technique requires that the surgeon be experienced in the placement of the pedicle screws in vertebrae. Before attempting this technique, surgeons are advised to attend a training course with a surgeon already experienced with the use of the device.

1 | DESCRIPTION

The Erisma-LP system is composed of rods fixed on the spine with pedicle screws. Moreover, this instrumentation includes transverse link which connects two rods altogether. The implants used in the Erisma-LP system are made of titanium alloy ISO 5832-3 or ASTM F136 as indicated by the symbol . It is essential to insert implants with instrumentation specifically designed for

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this purpose. For more description of the instrumentation it is necessary to read the technical documentation associated to the Erisma-LP product. The Erisma-LP implants must be assembled with new Erisma-LP components defined as being compatible with one another.

2 | INDICATIONS

General criteria and principles related to instrumented spinal surgery are applied here :

- Degenerative pathologies
- Spinal deformity
- Traumatology
- Tumor pathology

3 | CONTRA-INDICATIONS

- Pregnancy
- Infection
- Recognized allergies to titanium or titanium alloys

4 | SECONDARY AND POSSIBLE SIDE EFFECTS

- Pseudarthrosis
- Implant disassembly or implant failure
- Infection
- Allergy to materials used

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5 | CAUTIONS OF USE

Never reuse an implant, even in perfect state. Any implant which has been used, twisted, bent, implanted and then removed, even if appears intact, must be discarded. Use new implants routinely. This is indicated on the labeling of the implant package by the next symbol . Do not mix dissimilar materials.

- **Preoperative** : The surgeon must be fully conversant with all aspects of the surgical technique and know the indications and contra-indications of this type of implant. The surgeon must have acquainted himself before the operation with the specific technique for insertion of the product which is available from the manufacturer. As part of the pre-operative examination, the surgeon must check that no biological, biomechanical or other factors will affect the correct conduct of the operation and the post-operative period. An appropriate range of sizes must be available at the time of the operation.
- **Intraoperative** : The correct selection of the type of size of implant appropriate to the patient and the positioning of the implant are extremely important.

- **Postoperative :** Patients must be informed of the precautions to be taken in their everyday life to guarantee a maximum implant service life. It is recommended that a regular postoperative follow-up be undertaken to detect early signs of failure of the implants and to consider the action to be taken. Deterioration of the device after bone consolidation cannot be considered to constitute a dysfunction or degradation in the characteristics of the implant. The implant can be removed after the consolidation of the bone graft. A suitable rehabilitation program must be designed and implemented.

6 | HANDLING AND STORAGE

The handling of the Erisma-LP material must be done as seldom as possible and always with the

utmost care. The Erisma-LP implants (in their original packaging) must be stored with care in a clean and dry place. Do not expose the Erisma-LP implants to radiations or extreme temperatures. Should these requirements not be followed, reduced mechanical properties may occur which could lead to implant failure in some cases. Proper function of the surgical instruments specific to the Erisma-LP instrumentation should be verified prior to every surgical procedure.

7 | CLEANING / DECONTAMINATION / STERILISATION

The operations of cleaning, decontamination and sterilization must be realized for the implants and the instruments before and after use. All packaging must be removed before the operations of pre treatment, cleaning and

sterilization which will be implemented according to the recommendations good practices in force.

The implants must be sterilized before use. It is advisable to sterilize the implants in a steam autoclave according to the method used in hospitals and clinics (recommended values below).

Sterilization method	Autoclave
Cycle	Porous load (prion)
Temperature	134°C
Exposure time	18 Min (Minimum)

However we advise the users not following the recommended method to validate their methods by means of appropriate laboratory technique.

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