

Xtract-All[®] SPINE^{V.2}

Universal Spinal Implant Removal System

...Because Every Minute Counts[®]

Surgical Technique Guide




SHUKLA MEDICAL[®]

Universal Orthopedic Extraction Technologies

Contents

Introduction	3
List of Components	4-9
Cervical System	4
Thoracolumbar System	6
Instrumentation for Total Spine	8
Surgical Technique	10-13
Total Spine Techniques	10
Helicopter Method	12
Tips & Pearls for Cervical & Thoracolumbar Systems	13
Notes	14



Helicopter Socket with Screw Construct; page 12



Introduction

Intended Use

The Xtract-All® Universal Spinal Implant Removal System (S9SPINE) is designed to simplify spinal revisions.

Indications for Use

The Xtract-All® Spine (S9SPINE) is appropriate for any spinal revision case. For use by, or as directed by, a surgeon during spinal revision surgeries. The system includes over 130 implant drivers in a wide range of configurations and sizes, covering both standard and proprietary implant configurations. They can be used to remove hardware from virtually any spinal implant system.

Contraindications

The Xtract-All® Spine system is designed to be used when removing spinal hardware with intact screw heads. The system cannot be used with broken or stripped screws. For broken and stripped screw removal, please contact Shukla Medical Customer Service for information about the Xtract-All® Screw Universal Broken & Stripped Screw Extraction System (S9SCRW).

Warning and Precautions

Preoperative

- Clear x-rays and surgical notes may be used to identify manufacturer, brand, location, & condition of implanted hardware.
- The surgeon should be familiar with general principles of revision surgery and techniques for removal of implants.
- The instruments should be inspected for visible damage prior to use. Do not use the product if damage is suspected.
- Only validated cleaning and sterilization procedures should be used

Operative




- Proper handling and storage of the instrumentation is mandatory. Damage to the instrumentation may produce stresses and cause defects, which could become a focal point for failure.
- The surgeon should be cautious with spinal position change and/or excessive force exertion while removing implants using the instrumentation provided in the tray.
- All instrumentation has physical limits. Excessive force may result in instrument failure. It is recommended to maintain access to the Xtract-All® Screw Universal Broken & Stripped Screw Extraction System (S9SCRW) in the event that instrumentation fails.

Cleaning and Sterilization

For detailed sterilization instructions, please refer to industry standards ANSI/AAMI ST79:2012 & A1:2012 and ANSI/AAMI ST8:2001.

NOTE: All Shukla Medical surgical instruments require manual cleaning with a neutral pH cleanser. Open and disassemble all instruments, making sure to remove all contamination during cleaning. Instruments must be reassembled prior to sterilization. Machine washing is not recommended. Maintenance and care using an autoclaveable lubricant on movable parts is required to preserve the life of the instrument. For more cleaning, inspection, maintenance, and care tips, contact Shukla Medical directly.

For detailed cleaning and sterilization instructions, please visit www.ShuklaMedical.com/Sterilization.

EC	REP	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands	REF	S9SPINE	 Shukla Medical 8300 Sheen Drive St. Petersburg, FL 33709 USA	 CONSULT INSTRUCTIONS FOR USE	 NON-STERILE PRODUCT
----	-----	---	-----	---------	---	--	--



List of Components

Standard Drivers for Cervical Spine

Cervical Spine System

SCS034	Case 1: Cervical Spine System
SCS037	Lid: Cervical Spine System



Key

Case Location ID	Part #	Part Description
------------------	--------	------------------



Cruciform

49	SDR749	Cruciform 2.0mm
50	SDR750	Cruciform 2.5mm
51	SDR751	Cruciform 3.0mm
52	SDR752	Cruciform 3.5mm



Torx

56	SDR756	Torx T6
57	SDR767	Torx T7
58	SDR757	Torx T8
59	SDR768	Torx T9
60	SDR758	Torx T10
61	SDR759	Torx T15
62	SDR760	Torx T20
63	SDR761	Torx T25



Hex

1	SDR701	Hex 2.0mm
2	SDR702	Hex 2.25mm
3	SDR703	Hex 2.5mm
4	SDR704	Hex 2.75mm
5	SDR705	Hex 3.0mm
7	SDR707	Hex 3.5mm
8	SDR708	Hex 4.0mm
16	SDR716	Hex 7/64"
17	SDR717	Hex 1/8"
18	SDR718	Hex 5/32"
19	SDR719	Hex 3/16"



Flat

80	SDR780	Flat Small
81	SDR781	Flat Medium
82	SDR782	Flat Large



Phillips

83	SDR783	Phillips
----	--------	----------



Hex Sockets

24	SDR724	Hex Socket 3.5mm
25	SDR725	Hex Socket 4.0mm
37	SDR737	Hex Socket 1/8"
38	SDR738	Hex Socket 5/32"
39	SDR739	Hex Socket 3/16"
40	SDR740	Hex Socket 7/32"
46	SDR746	Hex Socket 1/2"



Saddle

84	SDR784	Saddle 4.0mm
----	--------	--------------



87	SDR787	Blade 4.0mm
----	--------	-------------



Hexalobe

145	SDR875	Hexalobe X10
146	SDR876	Hexalobe X15
147	SDR877	Hexalobe X20



Parts not shown to scale

Drivers for Proprietary Configurations



Aesculap®

95 SDR800 Aesculap® 5-Star



Alphatec Spine®

96 SDR801 Alphatec® Unlock Tool



Blackstone®

97 SDR802 Blackstone® 2-Prong



98 SDR803 Blackstone® Tri-Lobe



Biomet®

77 SDR777 Biomet® Pentalobe S15



Corin®

99 SDR804 Corin® Cervive 3-Prong



EBI®

104 SDR809 EBI® SpineLink® ACS



Interpore Cross®

76 SDR776 Interpore Cross® Pentalobe



101 SDR806 Interpore Cross® Unlocking Tool



Medtronic®

102 SDR807 Medtronic® Quad



Orthofix®

69 SDR769 Orthofix® Square 2.0mm

70 SDR770 Orthofix® Square 2.5mm

71 SDR771 Orthofix® Square 3.0mm

72 SDR772 Orthofix® Square 3.5mm



Stryker®

103 SDR808 Stryker Spine® 4-Prong



Zimmer®

105 SDR810 Zimmer Spine® Cervi-Lok®



107 SDR812 Zimmer® Nex-Link®



For the most up-to-date information on supported implants, use the X-Ray Search Engine™ Implant Compatibility Database.

The X-Ray Search Engine™ makes cross-compatibility identification easy. Enter keywords for any implant, manufacturer, or brand to find a list of all compatible Xtract-All® tools and related surgical techniques, x-rays, and videos. Our database is updated frequently as we are made aware of other implant systems.

Contact customer service or your sales representative to access the X-Ray Search Engine™.



List of Components

Standard Drivers for TL Spine

TL Spine System

SCS035	Case 1: TL Spine System
SCS038	Lid: TL Spine System



Key

Case Location ID	Part #	Part Description
------------------	--------	------------------



Cruciform

51	SDR751	Cruciform 3.0mm
52	SDR752	Cruciform 3.5mm
53	SDR753	Cruciform 4.5mm
54	SDR754	Cruciform 6.0mm



Torx

61	SDR759	Torx T15
62	SDR760	Torx T20
63	SDR761	Torx T25
64	SDR762	Torx T27
65	SDR763	Torx T30
66	SDR764	Torx T40
67	SDR765	Torx T45
68	SDR766	Torx T50



Hex

5	SDR705	Hex 3.0mm
7	SDR707	Hex 3.5mm
8	SDR708	Hex 4.0mm
9	SDR709	Hex 4.5mm
10	SDR710	Hex 4.7mm
11	SDR711	Hex 5.0mm
12	SDR712	Hex 6.0mm
13	SDR713	Hex 7.0mm
20	SDR720	Hex 7/32"
21	SDR721	Hex 1/4"



Flat

82	SDR782	Flat Large
----	--------	------------



Phillips

83	SDR783	Phillips
----	--------	----------



Saddle

84	SDR784	Saddle 4.0mm
85	SDR785	Saddle 5.0mm
86	SDR786	Saddle 6.0mm
87	SDR787	Blade 4.0mm
88	SDR788	Blade 5.0mm
89	SDR789	Blade 6.0mm



Pedicle Screws

138	SDR868	4-Prong Pedicle Screw
-----	--------	-----------------------



Hexalobe

146	SDR876	Hexalobe X15
147	SDR877	Hexalobe X20
148	SDR878	Hexalobe X25



Hex Sockets

26	SDR726	Hex Socket 5.0mm
27	SDR727	Hex Socket 5.5mm
28	SDR728	Hex Socket 6.0mm
29	SDR729	Hex Socket 7.0mm
30	SDR730	Hex Socket 7.5mm
31	SDR731	Hex Socket 8.0mm
32	SDR732	Hex Socket 9.0mm
33	SDR733	Hex Socket 10.0mm
34	SDR734	Hex Socket 11.0mm
41	SDR741	Hex Socket 1/4"
42	SDR742	Hex Socket 9/32"
43	SDR743	Hex Socket 5/16"
44	SDR744	Hex Socket 3/8"
45	SDR745	Hex Socket 7/16"

Parts not shown to scale



Drivers for Proprietary Configurations



Advanced Spine®

120 SDR850 Advanced Spine® Sleeve Nut



121 SDR851 Advanced Spine® Anchor



Aesculap®

124 SDR854 Aesculap® 4-Prong



Biomet®

77 SDR777 Biomet® Pentalobe S15



139 SDR869 Biomet® Polaris® 5.5



Corin®

128 SDR858 Corin® 3-Lok



Depuy®

100 SDR805 Depuy Spine® Outer Nut



129 SDR859 Depuy® Moss-Miami® Socket



130 SDR860 Depuy® Moss-Miami® Interrupted Socket



EBI®

132 SDR862 EBI® Square



133 SDR863 EBI® SpineLink® II Socket



134 SDR864 EBI® SpineLink® II



Interpore Cross®

125 SDR855 Interpore Cross® Hex 13.8mm



126 SDR856 Interpore Cross® Cap Nut Remover



127 SDR857 Interpore Cross® 2-Prong



76 SDR776 Interpore Cross® Pentalobe



Orthofix®

73 SDR773 Orthofix® Square 4.0mm



Theken Spine®

142 SDR872 Theken Spine® Driver



Spine Tech®

123 SDR853 Spine Tech® Octagon



Stryker®

135 SDR865 Stryker® Xia® Pedicle Screw



136 SDR866 Stryker® Xia® III Hexalobe



Synthes®

122 SDR852 Synthes® 11mm 12pt Hex



131 SDR861 Synthes® Click'X®



Zimmer®

106 SDR811 Zimmer® Incompass®



137 SDR867 Zimmer® Silhouette® Heptagon 1/2"



143 SDR873 Zimmer® Dynesys® L.I.S.



149 SDR879 Zimmer® Pathfinder® Hexalobe E7



US Spine®

144 SDR874 US Spine® 3-Prong




List of Components

Instrumentation for Total Spine

Key

Part #	Part Description
--------	------------------


Instrument Case



SCS036	Case 2: Spine Instrumentation
SCS039	Lid: Spine Instrumentation

Handles

MHRS0311	Ratcheting In-Line Handle
HD233-M01	Ratcheting T-Handle



Helicopter Sockets

SDR813	Helicopter Socket Small
SDR814	Helicopter Socket Medium
SDR815	Helicopter Socket Large



Extensions

SXN007	Extension Shaft 4"
SXN008	Extension Shaft 6"
SXN009	Extension Shaft 8"
SXN010	Extension Shaft 12"



Parts not shown to scale

Wrenches

- SWR003 Double Ended Wrench 7/32" x 9/32"
- SWR004 Double Ended Wrench 1/4" x 3/8"
- SWR005 Double Ended Wrench 5mm x 7mm
- SWR006 Double Ended Wrench 6mm x 10mm




Breaker Bar

- HD239-M01 Breaker Bar



Power Adapter

- SAD005 Power Adapter




Locking Pliers

- SWR008 Long Nose Locking Pliers 9"



Rod Grippers

- SWR009 Rod Grippers



Surgical Technique for Total Spine

1 Identification & Selection

Identify the spinal implant system from the surgical notes and X-rays. Select the appropriate Xtract-All® drivers. For the most up-to-date information on supported systems, contact customer service or your local rep for access to our X-Ray Search Engine™ Implant Compatibility Database.



- If a range of drivers is recommended or the spinal implant system cannot be identified, inspect the locking nuts & screws to visually select the most appropriate driver.
- If the appropriate driver cannot be identified, or a nut or screw is unable to be removed from the construct, the Helicopter Method may be used.

Note: For broken and stripped screws, please contact Shukla Medical Customer Service for information about the Xtract-All® Screw (S9SCRW) Universal Broken & Stripped Screw Extraction System.

Pre-op Planning: For assistance identifying implants and determining compatible drivers, please contact Shukla Medical Customer Service and let our team of experts help you.

2 Assemble Driver

Insert the selected driver into the appropriate Extension Shaft (SXN007-SXN010) if needed. Insert Extension Shaft into Ratcheting Screwdriver Handle (HD233-M01 or MRHS0311). Rotate handle to change between ratcheting mode.

- Center:** Ratcheting disabled
- Forward:** Insert screws
- Back:** Remove screws



Multiple handle styles are available depending upon surgeon preference and desired level of torsion.

Ratcheting In-Line Handle

MRHS0311

Suitable for most implant removals. Ratcheting mode switches between forward, backward, or neutral/locked.



Ratcheting T-Handle

HD233

Suitable when additional torque is required during manual implant removal. Ratcheting mode switches between forward, back, or neutral/locked. Use while locked for the Helicopter Method (facing page).



Breaker Bar

HD239-M01

A breaker bar is included in case of difficulty due to well-fixed screws. Using the breaker bar can generate significant torsional force that may not be optimal in some spinal procedures. Use with caution.



3 Rod Removal

- a. Remove locking nuts with assembled screwdriver (Fig. A, panels 1 & 2).
- b. Stabilize & remove spinal rods using Rod Gripper (SWR009) & Long Nose Locking Pliers (SWR008) (Fig. A, panels 3 & 4).




Figure A

4 Screw Removal

Select the appropriate driver. Assemble screwdriver as per Step 2.

- **Note:** For screws that do not have an internal configuration at the bottom of their uni-axial screw, use the blade or saddle drivers (SDR784-SDR815). If available drivers do not fit, reassemble the locking nut & proceed to use the *Helicopter Method* (pg 12).

Remove screw using screwdriver assembly. Locking Pliers (SWR008) may be used to aid with removal. (Fig. B)



Figure B

OPTIONAL: *The Helicopter Method*

The Helicopter Method facilitates total screw construct removal by rotating the pedicle screw while still attached to the rod using a Helicopter Socket (SDR813-SDR815).


The screw construct consists of the screw, rod, & locking cap.

Indications for the *Helicopter Method*:

- If the correct driver cannot be identified
- If any cap, nut, or screw is fixed so tightly that it cannot be removed
- If the locking cap or screw is stripped

1 Cut rod on either side of tulip

- Approximately 5mm of rod should remain extending from sides of tulip head.
- Rod cut length must be long enough to engage with helicopter socket, but short enough to minimize damage to surrounding live tissue as screw construct rotates.



2 Assemble driver with Helicopter Socket

- Select Helicopter Socket (SDR813-SDR815) that best fits over tulip head.
- Connect socket to an Extension Shaft (SXN007-SXN010), then attach socket assembly to T-Handle (HD233-M01).
 - Ratcheting mode must be locked or set to reverse.



3 Use Helicopter Socket to remove screw construct

- Screw construct must be fully assembled in order for the *Helicopter Method* to be effective.
- Place Helicopter Socket over tulip and rod, so that rod is engaged in socket grooves.
- Turn counterclockwise until screw construct backs out.
 - If additional torque is needed, attach Breaker Bar (HD239-M01) to extension and turn counterclockwise.



Tips & Pearls

Cervical System


The Cervical Case includes 56 drivers covering standard and proprietary implant configurations.



Plate Removal via the *Levitating Plate Method*

Bone screws secured with locking o-rings

If screws holding plate are secured with locking o-rings, plate may be removed by loosening all bone screws incrementally so that they lift the plate evenly from the surface of the bone all together.



Thoracolumbar System

The TL Case includes 76 drivers covering standard and proprietary configurations.



Rod Removal

Well-fixed locking nuts

If additional torque is required, use the breaker bar (HD239-M01). Using the breaker bar can generate significant torsional force that may not be optimal in some spinal procedures. Use with caution.



For more information about the
Xtract-All® Spine (S9SPINE) Universal Spinal Implant Removal System
Call us at 888-4-SHUKLA (888-474-8552) or visit
www.ShuklaMedical.com



...Because Every Minute Counts®