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Characteristic

The Occipitocervical Spine Fixation System is a set of implantable metallic components (medical product), complementary and compatible, developed for the application in spine surgery, with the purpose of promoting stabilization, or of the occipitocervical region and the cervical segment, according to the indication for use.

Material

Titanium alloy according to ASTM F-136

This surgical technique was developed by:

Technology Division - Baumer S.A.

Indications

The indications of the Occipitocervical Spine Fixation System are for cases of correction, stabilization and fixation in the occipital / cervical region and in the posterior cervical segment of the spine, due to diseases, such as:

- Traumatology (including fractures);
- Deformities;
- Degenerative and traumatic pathologies;
- Tumors;
- Sequelae of inadequate or unsuccessful treatments for the pathologies described above.

Contraindications

The contraindications below must be strongly considered by the responsible surgeon. Such contraindications include, but do not limit the use of the Spine Fixation System - Occipitocervical Fixation, in cases of

- Recent history of systemic or localized infection;
- Signs of local inflammation;
- Bone immaturity;
- Immunosuppressive pathologies;
- Mental disease;
- Alcoholism or drug abuse;
- Metastatic tumors spread widely in the adjacent vertebral bodies;
- Pregnancy;

- Severe osteopenia;
- Severe instability due to the lack of structural tissue integrity;
- Severe osteoporosis;
- Morbid obesity;
- Fever or leukocytosis;
- Sensitivity or allergy to the metal of the medical product;
- All those known for the use of column fixation systems;
- Patients unwilling or unable to follow the post-operative care and instructions;
- Any medical or surgical condition that could compromise the success of the surgical procedure;
- Every patient should be alerted to the risks of walking without the aid of a crutch or cane, before biological fixation occurs, complete recovery and adaptation to the implantable metallic components used.

Introduction

The Occipitocervical Spine Fixation System is composed of a series of implants, such as: 3 and 4-hole Occipitocervical Bar, Longitudinal and Cross Bar, Articulated Screw (Poliaxial) with Total Thread and Partial Thread, Spongy Screw, Screw Support, Long and Short Angled Hook, Occipital Screw and Nut, Upper Nut, Longitudinal and Square Connectors, Offset Connectors, among others, whose main purpose is to align, stabilize, correct and fix the segments of the occipital / cervical and posterior cervical spine vertebral.

The Articulated Screws (Poliaxial) with Total and Partial Thread, have a thread diameter of Ø3.5mm and Ø4.0mm and guarantee good anchorage close to the vertebral bone tissue. The polyaxial head of the articulated screw allows placement in angled positions, providing greater freedom to the surgeon in choosing the best position for the screw. All screws are self-cutting, eliminating the use of the tap.

Occipitocervical bars in versions with 3 and 4 holes allow the doctor to choose the model that best adapts to the patient's profile and the need for fixation. The holes are oval, offering freedom in the placement of the screws and a fine adjustment of the final assembly.

Longitudinal bars are available in diameter Ø3.5mm, with lengths from 30mm to 260mm.

The Crossbars can be fixed through connectors, after the final assembly of the set and provide stability and rigidity to the system.

- . Angled Hooks are also available in the Short and Long models, right and left side, for fixing and anchoring close to the atlas (C1).
- . The longitudinal and square connectors allow the connection of the occipital assembly to the cervical segment, being able to extend the assembly to the thoracic segment of the spine, with the proper complement of the Pedicular Fixation System.
- . The closed OC Offset Connectors make it possible to join an articulated screw or hook to the longitudinal bar. This connection occurs in a displaced manner, in order to avoid excessive molding of the longitudinal bar. They are 10, 15 and 20mm long.
- . Fixation in the occipital region can be done through Reverse occipital screw / Occipital nut or through Ø3.5mm and Ø4.0mm spongy screws.
- . The instruments for handling and inserting the implants are simple, compact and allow a safe placement of the components.

Technique


Implant Selection

It is essential that preoperative planning is performed to avoid unexpected occurrences during the surgical procedure. This procedure should only be performed by surgeons qualified and trained with Baumer instruments.

ATTENTION: Only use instruments with the Baumer logo. The company has a periodic instrument maintenance program to ensure that distributors provide good quality instruments. If you receive instruments with problems, contact SAC - Customer Service, by e-mail baumer@baumer.com.br, informing the code and lot engraved on the part, as well as the name of the distributor.

1. Access route

Using the standard surgical approach, expose the vertebral region to be operated.


2. Initial drilling

Determine the entry point and the screw's path through the initiator IN7.120.001 create a pilot hole where you want to insert the screw.


After the initial drilling, use the Adjustable Drill Guide CH.75.C and the drill 900.3008 to drill the vertebral body corresponding to the length of the screw to be implanted. If necessary, use the pedicle probe IN.5 to check the position of the cavity.

Drill in all the vertebral bodies where you want to place the screws.

3. Depth measurement


After drilling with the drill, check that the hole depth is adequate for the placement of the respective screw. Use the ME.17 depth gauge to make this measurement.

4. Introduction of the Articulated Screw

Fit the predetermined articulated screw to the tip of the IN6.041.003 Screwdriver, then turn the lock clockwise to thread the screw and fix it to the Screwdriver.

Insert the screws into the holes made by the drill. To remove the Screwdriver that is attached to the screw, simply turn the locking device counterclockwise. Detail of hinged screws inserted.


5. Placement of the Support and Screw for Lateral Mass

Mount the screw on the IN6.041.004 introducer wrench and secure it with the clamp. Then mount the screw support 860.0.D / E, 860.15.D / E or 860.35.D / E with an angle of 0°, 15° and 35° on the IN7.600.013 threaded wrench. Insert the IN7.600.012 guide Screwdriver into the slot in the support and give the final tightening to the set. Then position the support on the vertebra (lateral mass), and insert the screw.

Detail of hinged screws placed and screw holders for lateral mass.

6. Molding of the Longitudinal Bar

By inserting the templates TT.32.8 / 16 / 25 in the articulated screws and implanted screw support, it is possible to obtain a pre-molding of the longitudinal bar. The molding of the 863.REF longitudinal bar must be made with the IP.66 modeler or with the IP.65 modeler, using the template that is already molded as a guide.

7. Placement of the Longitudinal Bar

After the longitudinal bar has been previously molded, place it on the articulated screws and screw supports with the aid of the IN7.590.002 bar holder clamp.

! Attention: The bar must be fully accommodated. If necessary, use the GR.128 bar pusher to assist in the procedure.

8. Placement of the Lock Screw

Mount the 864 locking screw on the IN6.041.004 wrench. To facilitate the placement of the locking screw on the hinged screws and screw support, the screwdriver must be turned counterclockwise until you realize that the locking screw has engaged the thread. There is a "CLICK" sound, and then turn the screwdriver clockwise to start the installation. In this step, it is not necessary to give the final tightening, just ensure that the locking screw has been placed correctly (check the alignment) and that the longitudinal bar is properly fixed. Detail of the locking screws placed.

! The correct alignment of the locking screw must be ensured before starting the fitting.

Occipitocervical fixation

To ensure a safe procedure, it is recommended that the thickness of the occipital bone be measured with the help of x-rays or other imaging possibilities. The thickness will indicate the size of the occipital screws that will be implanted later.

9. Bar Molding

By inserting the TT.32.8/16/25 templates in the articulated screws and implanted screw support, it is possible to obtain a pre-molding of the bar. The molding of the occipitocervical bar must be done with the IP.66 modeler or with the IP.65 modeler, using the template that is already molded as a guide.

Attention: Only LONGITUDINAL and OCCIPITOCERVICAL BARS can be mechanically shaped or molded, and the limit for the molding of the bars is the normal curvature of the column. It is worth noting, however, that the bars cannot be shaped / bent at sharp angles, scratched or deformed. In addition, the molding can never be reversed, that is, once molded, the health product cannot be molded back to its original shape, which can cause fracture and failure of the fixation system.

10. Placement of the Occipitocervical Bar / Marking

After the occipitocervical bar has been previously molded in the articulated screws and screw supports with the aid of the IN7.590.002 bar holder forceps, make sure that the bar is well supported and aligned in the occipital region. Make a mark on the bone where the reverse occipital screw will be placed, remove the bar and complete the hole / slot mark. This marking will serve as a guide for making the final hole / tear in the occipital region.

11. Drilling for Occipital Fixation


After marking, drill the hole that will allow the 865.REF reverse occipital screw head passage.

NOTE: Fixing the bar / plate to the occipital can also be done with the 3.5mm 860.REF screw. Drill the initial hole using the adjustable drill guide CH.75.C and drill 900.3008, resting on the hole in the occipitocervical bar. Then insert the 3.5mm screw to lock. This fixation is less aggressive and more than 02 screws can be placed in each occipitocervical bar.


12. Placement of the Reverse Occipital Screw

Select the 865.REF reverse occipital screw suitable for the bone structure. Thread the screw on the CH.101 guide rod. Then, insert the screw head into the occipital hole.

13. Placement of the Occipitocervical Bar on the Reverse Occipital Screw

In this step, the reverse occipital screw remains threaded on the guide rod and positioned in the occipital. Place the occipitocervical bar on the guide rod. Then, mount the lock nut 865 on the CH.101.C. Place the wrench with the nut on the guide rod that is threaded on the screw, and start tightening the assembly consisting of: reverse occipital screw / occipitocervical bar / nut.

14. Final Tightening of the Occipitocervical Bar on the Reverse Occipital Screw

Remove the introducer wrench and the guide rod that is threaded on the reverse occipital screw. Mount the final torque wrench CH.101.B over the anti-torque wrench CH.101.A. Insert the anti-torque wrench into the hole and slot of the reverse occipital screw to prevent the occipital screw from turning during the final tightening. Give the final tightening by turning the tightening screwdriver clockwise until the assembly is fixed.

Detail of the occipitocervical bar placed in the occipital region and in the posterior cervical region.

To place the bar on the other side, follow the same sequence used previously.

15. Placement of the Angled Hook - Complementary Fixation

Some assemblies require additional fixation. In these cases, you can use the short hooks or long hooks for the right and left side.

After selecting the hook 862.A/B/C/D, mount it on the IN7.600.013 screwdriver. Fit the IN7.600.012 guide wrench into the slot of the angled hook and tighten the set. Then place the angled hook fixing it on the posterior arch of the C1 vertebra, and then place the lock nut to fix the angled hook on the occipitocervical bar.

Final Fixation

16. Distraction and Compression System

Distraction and compression movements are possible with the forceps CH.103 and CH.103.A.

NOTE: To obtain these movements, one of the articulated screws, angled hooks or screw support must be tight while the other is only placed, but without tightening. After distraction or compression,

tighten the locking screw to fix the set.

! It is worth noting that excessive tightening should be avoided. The final tightening will be done, later, with the aid of a Torque Wrench and an Anti-Torque Wrench

17. Final tightening

After aligning the vertebra, perfectly accommodating the bar on the base of the tulip and ensuring the desired correction, the final tightening of the locking screws must be performed.

For this, it is necessary to position the Anti-Torque Screwdriver IN6.001.003 next to the screw tulip (upper / proximal portion of the screw). Mount the final tightening wrench IN6.011.001 on the IN7.600.014 torque wrench, insert it into the Anti-Torque Wrench, and fit the hex end of the wrench to the hex of the locking screw. Torque SLOWLY until you hear a “click”. Do this for all screws / hooks.

ATTENTION: The screws of the Occipitocervical Fixation System are designed for a maximum final tightening of 2.0 Nm. The tightening movement must be STOPPED when a click is heard (which indicates that the proper torque has been reached), so that no breakage or damage occurs. of the torque wrench or screw deformation. It is not necessary to apply more torque than specified, to “guarantee” the final tightening, because instead of increasing the performance of the system, the opposite may occur: deformation of the implant components, loosening of the locking screw or other type of failure in the system due to excess torque

ATTENTION: The torque wrench must never be used to make the initial tightening and to remove the locking screws. For this task, the instrument kit provides the physician with another 600.1000 "T" 600.1000 that can be used with the IN6.011.001 wrench or the other wrench already informed in the surgical technique.

The Anti-Torque Wrench MUST be used during the final tightening. Failure to observe this step may compromise the surgical procedure and its results. The Anti-Torque Wrench has two important functions: 1. It allows the torque wrench to align with the clamping axis; 2. It allows the optimization of the torque required to block the implant assembly. ATTENTION: If the Anti-Torque Wrench is not easily removed from the tulip, the longitudinal bar may not be fully accommodated at the base of the screw tulip.


NOTE: According to the Instrument Maintenance Policy and aiming to ensure that the torque wrench (the most critical item of Occipto instruments) is always calibrated, the IN7.600.014 instrument must be sent for calibration and maintenance every 6 months in 6 months.

18. Placement of Cross-Links

Use the 866 bar / bar connector with the 866.REF crossbar to do the assembly.


Instruments - K.CG.5.001

Components	Description
CG.5.001	TRAY FOR OCCIPTOCERVICAL INSTRUMENTAL (STAINLESS)
IN7.600.012	HOOK AND SCREW BRACKET INTRODUCER
IN6.041.003	ARTICULATED SCREW INTRODUCER WRENCH
IN6.041.004	SCREW INTRODUCING WRENCH
IN7.120.001	PEDICLE INITIATOR
IN7.600.013	HOOK AND THREAD SCREW SUPPLY INTRODUCER
ME.17	CORTICAL METER
CH.101.B	WRENCH TO TIGHTEN OCCIPTAL NUT
CH.101	OCCIPTAL SCREW GUIDE
CH.101.C	WRENCH TO PLACE THE OCCIPTAL NUT
CH.75.C	ADJUSTABLE DRILL GUIDE FOR OCCIPTOCERVICAL
900.3008	DRILL Ø2.0MM FOR OCCIPTOCERVICAL
CH.101.A	WRENCH TO HOLD OCCIPTAL SCREW
GR.128	BAR PUSH FOR OCCIPTOCERVICAL
IN.5	PEDICLE PROBE FOR OCCIPTOCERVICAL

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