

THE LEADING CHOICE FOR
CERVICAL MOTION PRESERVATION¹

Mobi-C[®]

Cervical Disc



HIGHRIDGE

Motion Preservation

Mobi-C is an artificial disc, made to move with your patient as they move their neck, simulating motion similar to that of the natural cervical spine.

- **Superior to fusion for two-level disc replacement:**²

Mobi-C Cervical Disc is the first cervical disc in the United States approved to treat more than one level of the cervical spine.

- **One-step bone sparing insertion:** Mobi-C implantation requires no additional exposure or operative steps for screw or keel placement. No invasive drilling or chiseling are required.

- **Pre-assembled implants:** Mobi-C is delivered pre-assembled on a disposable PEEK cartridge, which fits easily on the implant inserter, saving operative steps. The PEEK Cartridge allows a radiolucent view of the implant for optimal positioning.




Over 225,000 Mobi-C Discs have been implanted in 25 countries since 2004.

Anatomically designed.
Clinically compelling.



WATCH
IT WORK



MOBILE BEARING TECHNOLOGY

The controlled mobility of the patented mobile core is the foundation of Mobi-C. With vertebrae and neck muscle movement, the Mobi-C implant is free to twist and slide left-to-right, and front-to-back, as well as rotate.

For more information, visit highridgemedical.com



HIGHRIDGE

CONFIDENCE THAT MOVES

1. Data on file. Based on available market data at the time of this publication.
2. 10-year data: Published 10-year outcomes of 1- and 2-level Mobi-C patients available.

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of Highridge Medical, LLC or an affiliate; and all products are manufactured by one or more of the spinal subsidiaries of Highridge Medical, LLC and marketed and distributed by its authorized marketing partners. Please refer to the Instructions for Use, package labeling, and applicable Surgical Technique Guide. Product clearance and availability may be limited to certain countries/regions. This material does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Highridge Medical. HM0221 REV A 05/25 ©2025 Highridge Medical, LLC.