


TALON™ TLIF CAGE

► Indications



- Low grade I or II Spondylolisthesis
- Severe or advanced degenerative disc disease
- Instability of the anterior column in association with posterior pathology


► Specific Advantages

- TALON™ TLIF cage makes to change the angle of the implant and to be inserted easily
- TALON™ TLIF cage is available in 4° lordosis
- Minimal implant width and tapered tip to facilitate the insertion through a small opening
- Anatomic shape to conform to the anterior vertebral body
- Bullet shape



IMPLANT DESCRIPTION

ARTICLE .NO	H(mm)	W(mm)	L(mm)	A
TT0826	8	11	26	2°
TT0926	9	11	26	2°
TT1026	10	11	26	2°
TT1126	11	11	26	2°
TT1226	12	11	26	2°
TT1326	13	11	26	2°
TT1426	14	11	26	2°
TT1526	15	11	26	2°
TT0832	8	11	32	2°
TT0932	9	11	32	2°
TT1032	10	11	32	2°
TT1132	11	11	32	2°
TT1232	12	11	32	2°
TT1332	13	11	32	2°
TT1432	14	11	32	2°
TT1532	15	11	32	2°
TT0835	8	11	35	2°
TT0935	9	11	35	2°
TT1035	10	11	35	2°
TT1135	11	11	35	2°
TT1235	12	11	35	2°
TT1335	13	11	35	2°
TT1435	14	11	35	2°
TT1535	15	11	35	2°



INSTRUMENT VIEW

