



# ZAVATION

## ▶ Z-Link<sub>PC</sub> (PC-Posterior Cervical)

### ▶ Z-Link<sub>PC</sub> Benefits and Features

- Variety of dominoes gives surgeons options to connect rods without disturbing the construct already in place
- Screw head color matches drill & tap
- Same torque applied to every implant in PC reducing the risk of mistakes for surgeons
- Offers curved, straight, and transition rods giving the surgeon options



<b>Polyaxial Screws (25 lbs/in)</b>	<u>Diameter:</u>	<u>Length:</u>
Standard screws	3.5mm, 4.0mm, 4.5mm	10-52mm
Smooth shank screws	3.5mm, 4.0mm	20-50mm
<b>Rod Connectors (dominoes)</b>	<u>Sizes:</u>	
Parallel	• 3.5x3.5mm	
Inline	• 3.5x4.5mm	
Open	• 3.5x5.5-6.0mm	
<b>Rods</b>	<u>Sizes:</u>	
3.5mm Straight	20-120mm 10mm increment	150mm, 200mm, 240mm
Transition	60, 80, 100, 120, 175, 240mm (19mm length of 5.5mm diameter)	450mm (250mm of 5.5mm diameter)
Curved	20-120mm 10mm increments	Curvature: 24cm radius (9.5")
<b>Cross Links</b>	<u>Sizes:</u>	
28-68mm Rod Spacing (With a snap on attachment feature)	• 28-34mm • 34-45mm • 45-68mm	
<b>Hooks</b>	<u>Sizes:</u>	
Straight, left, & right	4.5mm and 6.0mm openings	
<b>Offsets</b>	<u>Lengths:</u>	
Open & Closed	10mm and 14mm	

220 Lakeland Parkway — Flowood, MS 39232 — 601.919.1119 — [info@zavation.com](mailto:info@zavation.com)

[WWW.ZAVATION.COM](http://WWW.ZAVATION.COM)

\*The information contained in this document is Zavation Medical Products, LLC. Proprietary Information. It shall not be used, disclosed to others, or reproduced without the expressed written consent of Zavation Medical Products, LLC. For full instructions for use, visit <http://zavation.com/general-information/>