


# ZYFIX™

## FACET FUSION SYSTEM



- ❑ Promotes stabilization of the facet
- ❑ True fixation system incorporating a titanium compression screw
- ❑ Unique fenestrations encourage fusion with use of autogenous bone graft
- ❑ Minimally invasive procedure
- ❑ Multi-level use



An Xtant Medical Company



## Dual-Thread Screws



## Single-Thread Screws




❑ 6.0mm diameter screws

❑ Simple Instrumentation

❑ Multiple surgical approaches

❑ Different thread pitches provide increased superior to inferior facet compression

\*Optional 11mm washer available



❑ Lower profile than traditional pedicle screws

❑ Anti-migration indent

❑ Top threads provide additional locking

❑ Fenestrated for autogenous bone graft

❑ Cannulated for percutaneous approach

### Zyfix™ Screws and Washer

Item Number	Description
X076-0020	6.0mm x 25mm Dual Thread Screw
X076-0021	6.0mm x 30mm Dual Thread Screw
X076-0022	6.0mm x 35mm Dual Thread Screw
X076-0023	6.0mm x 40mm Dual Thread Screw
X076-0024	6.0mm x 45mm Dual Thread Screw
X076-0029	6.0mm x 20mm Single Thread Screw
X076-0037	11mm Washer



### X-spine Systems, Inc.

452 Alexandersville Rd., Miamisburg, OH 45342  
 Phone: 800-903-0640 • Direct: 937-847-8400 • Fax: 937-847-8410  
[www.x-spine.com](http://www.x-spine.com)

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: X-spine products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

X-spine®, the X-spine logo®, and ZYFIX™ are trademarks or servicemarks of X-spine Systems, Inc.

Patents Pending.

© 2016 X-spine Systems, Inc., All rights reserved.

F-1000.28 Rev C 6/16

