



NILE™ Alternative Fixation



PRODUCT OFFERING

NILE™ ALTERNATIVE FIXATION CLAMPS



Part Number	Description
5401-84500	Ø4.5 mm Clamp
5401-85500	Ø5.5 mm Clamp
5401-84500H	Ø4.5 mm Rail & 6.0 mm Rod & 6.35 mm Rod Clamp
5401-85500H	Ø5.5 mm Rail Clamp


NILE™ ALTERNATIVE FIXATION STERILE PACKED BANDS





Part Number	Description
5416-F03730-SG	Narrow Band, 3 mm
5416-F04730-SG	Wide Band, 4 mm

NILE™ ALTERNATIVE FIXATION SPINAL SYSTEM

IMPLANTS DESIGNED TO PERFORM



- Four Clamp Options That Accommodate Six Rod Sizes/Types
- Audible Click & Tactile Feedback Upon Rod Engagement
- Smallest Run on the Rod in the Industry
- Band & Rod Locking Mechanisms are Independent of Each Other, Allowing Intra-Operative Flexibility




- Stainless Steel Leaders Provide Familiar Handling Characteristics
- Two Leader Length Options to Accommodate Various Techniques: 50 mm & 110 mm
- Wide & Narrow Band Options For Different Sizes of Anatomy
- Color-Coded to Indicate Orientation
- Smooth Surface Finish to Minimize Bone Surface Abrasion During Passage
- Biocompatible Woven Polyethylene Terephthalate (PET) Provides Strength & Maintains Structure

INSTRUMENTS DESIGNED FOR SIMPLICITY & CONTROL

- No Assembly Required
- One Tensioner Fits All Clamp Sizes

*Over 75% Smaller Than the Leading Competitor




- Adjustable travel distance allowing for large reduction
- Allows for controlled, sequential reduction

*Compared to *Universal Clamp*


CONSTRUCT TO PROVIDE STABILITY

STATIC PULL-OUT VS. UNIVERSAL CLAMP



The NILE Alternative Fixation System is smaller, but more robust & stable.

DYNAMIC COMPRESSION VS. UNIVERSAL CLAMP



Compared to NILE Alternative Fixation 5.5 mm Clamp and Wide 4 mm Band



NILE™
ALTERNATIVE FIXATION



www.K2M.com

©2015 K2M, Inc. All rights reserved.
K2-54-7004-01 Rev. 0
Actual device color may vary.
Consult product catalog for details.



K2M, Inc.
751 Miller Drive SE
Leesburg, Virginia 20175
USA

PH 866.K2M.4171 (866.526.4171)
FX 866.862.4144



Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands
PH +31.70.345.8570
FX +31.70.346.7299

