


IRENE

A Company of
Naton Medical Group

MIS PLUS

Thoracolumbar Spine Fixation System



OP-time less

Blood loss less

Fixation plus

 **NATON**

Features:

Minimally invasive surgery available smaller incisions, shorter recovery time

Round groove on both sides are used to connect the reduction forceps

Extension break-off tap length is 80cm, which protrudes from skin and simplifies matching instruments

Optional bone cement holes design may provide stronger fixation for the patients with osteoporosis

Increase the inner width of U-shaped extension tap, easy for rod insertion

Reverse angle thread design avoids tap break-off and thread damaging

Titanium Rod (MIS Plus)



Diameter (mm)	Length (mm)
5.5	50,55,60,65, 70,75,80,85, 90,95,100,110, 120,130,140,150

Polyaxial Reduction Screw



Diameter (mm)	Length (mm)
5	30,35,40,45
5.5	30,35,40,45
6	30,35,40,45
6.5	30,35,40,45
7	30,35,40,45

Indications:

- ▶ Thoracolumbar Spine Tumors
- ▶ Thoracolumbar Fractures
- ▶ Thoracolumbar Deformity
- ▶ Thoracolumbar Degenerative Disc Disease
- ▶ Thoracolumbar Revision Surgery