







COMPLETE SET OF SIZE FOR OPEN, MIS AND DEFORMITY SURGERY

BEST MECHANICAL PROPERT

COMPACT AND SOFISTICATED INSTRUMENTS SET

SCREW WITH SELF-TAPPING TIP FOR INSERTION EVEN WITHOUT TAPPING TOOL

ALWAYS WITH DOUBLE THREAD TO IMPROVE THE SEAL

SINGLE AND STERILE PACKAGING

ALWAYS CANNULATED FOR CEMENT AND BONE SUBSTITUTE

System for thoracolumbar vertebral stabilization

Misano is the first thoracolumbar stabilization system that involves the use of a single pedicle screw for the treatment of degenerative, traumatic or deformity pathologies.

This allows both open and percutaneous surgical techniques. The same screw can be cemented or convey bone substitutes.

The screw is self-tapping and can be introduced without the use of the tap, and has a specific double thread that guarantees rapid implantation and high pedicle sealing.

The system includes a range of hooks, dominoes and connectors as well as the same cannulas of the minimally invasive instruments that can be used for the treatment of deformities.







Main features








Instrument design

Clover has invested heavily in the design and care of the instruments with the aim of creating ergonomic, functional and compact instrumentation. Designed for the surgeon and the operating equipe.





MONOAXIAL SCREW CANNULATED - FENESTRATED	D5.5 D6.5 D7.5 D8.5	FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55
POLIAXIAL SCREW CANNULATED - FENESTRATED	D4.5 D5.5 D6.5 D7.5 D8.5	FROM L25 TO L40 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55
TITANIUM PRECURVED ROD	D5.5	FROM L25 TO L100
TITANIUM ROD	D5.5	FROM L110 TO L500
CO-CR ROD	D5.5	FROM L110 TO L500
CLAMPING CAP FOR SCREWS AND HOOKS		
DOMINOES		
LONGITUDINAL CONNECTION ELEMENT		
MODULAR CROSSLINK		S - M - L - XL
ANGLED OFFSET		FROM L20 TO L60
OFFSET		FROM L20 TO L60
LUMBAR WIDE HOOK		FROM SIZE 7 TO SIZE 11
PEDICLE HOOK		FROM SIZE 5 TO SIZE 9
TORACHIC LAMINAR HOOK		FROM SIZE 5 TO SIZE 9
OLBIQUE HOOK - DX/SX		
OFFSET HOOK - DX/SX		
LUMBAR NARROW HOOK		FROM SIZE 7 TO SIZE 11































Clover Orthopedics s.r.l.

Via Gadames n. 57/7, c.a.p. 20151 Milano

M. info@cloverorthopedics.com

W. cloverorthopedics.com

T. +39 02 457 902 31

F. +39 02 457 902 66



