

**Sürekli, yenilikçi ve
güvenilir fikirler**

*Continuous, innovative
and reliable ideas*



artrofiks
spinal

Shark Mobile Disc Prosthesis



www.artrofiks.com



Toothed surface for stable adherence to lower and upper vertebra

Code	Description	Size
CD12-04 ~ 16-09	Shark mobile disc prosthesis	12-04 ~ 16-9.0
CDP12-40 ~ 13-90	Shark mobile disc prosthesis Peek	12-0.4 ~ 13-9.0

Description

Artrofik's Mobile Disc Prosthesis are used for the purpose of keeping mobility after the discectomy fusion operation. It allows all motions in the directions of rotational and flexion extension. Interior of the implant are coated with a special alloy in order to keep continuous mobility and to reduce corrosion and friction quotient. Implant allows (+,-) 10° flexion/extension, (+,-) 10° lateral bending and (+,-) 8° rotation. Artrofik's Mobile Disc Prosthesis are in 2 forms as Disc prosthesis with Peek inner and Disc Prosthesis with Titanium inner.

Usage Area

- * MRI-CT compatible titanium material
- * Toothed nature that provides coherence to superior and inferior surfaces
- * Compatibility to Smith Robinson approach
- * Allowing insertion and stabilization with a single instrument
- * Various size option for different anatomies



Artrofik's Medikal San. ve Tic. Ltd. Şti.

Batı Sitesi Mah. Tahsin Kahraman Cad. No: 42

Yenimahalle Ankara TÜRKİYE

Tel : +90 312 219 77 51 - 219 69 91

Fax : +90 312 219 41 12