

# AURA ANTERIOR CERVICAL PLATE SYSTEM



**Perfect for  
cervical spine  
stabilization**

## **AURA Anterior Cervical Plate System HIGHLIGHTS**

- ☞ Titanium raw material
- ☞ Cervical Plate system for Anterior stabilization
- ☞ Can be used for one or more segments
- ☞ Screw locked
- ☞ Quite thin plate
- ☞ No height on the plate



# AURA ANTERIOR CERVICAL PLATE SYSTEM




## PURPOSE OF USAGE

AURA Anterior Cervical Plate System is used in anterior cervical internal fixation procedures to ensure cervical spine stabilization in the event of tumor, trauma deformity, degenerative disc disease and other forms of cervical instability in the spine. Fixation of the AURA Anterior Cervical Plate System is achieved by inserting bone screws into the cervical spine using the anterior approach.



## INDICATIONS

- ♀ Degenerative diseases of the cervical spine
- ♀ Infectious or inflammatory pathologies
- ♀ Trauma
- ♀ Cases undergoing corpectomy
- ♀ Degenerative disc diseases
- ♀ Tumor
- ♀ Pseudarthrosis





# AURA ANTERIOR CERVICAL PLATE SYSTEM

## AURA PLATE (4 HOLE)

Size	Reference Number
22mm	ACP02204
24mm	ACP02404
26mm	ACP02604
28mm	ACP02804
30mm	ACP03004
32mm	ACP03204
34mm	ACP03404
36mm	ACP03604



## AURA PLATE (6 HOLE)

Size	Reference Number
38mm	ACP03806
40mm	ACP04006
42mm	ACP04206
44mm	ACP04406
46mm	ACP04606
48mm	ACP04806
50mm	ACP05006
52mm	ACP05206
54mm	ACP05406



## AURA PLATE (8 HOLE)

Size	Reference Number
56mm	ACP05608
59mm	ACP05908
62mm	ACP06208
65mm	ACP06508
68mm	ACP06808
71mm	AU1071
72mm	ACP07208





# AURA ANTERIOR CERVICAL PLATE SYSTEM


## AURA PLATE (10 HOLE)

Size	Reference Number
68mm	AU1168
72mm	AU1072
76mm	ACP07610
80mm	ACP08010
84mm	ACP08410
88mm	ACP08810
92mm	ACP09210
96mm	ACP09610
100mm	ACP10010




## AURA SCREW

Size	Reference Number
Ø4.0X12	AUR4012
Ø4.0X14	AUR4014
Ø4.0X16	AUR4016
Ø4.0X18	AUR4018
Ø4.0X20	AUR4020
Ø4.5X12	AUR4512
Ø4.5X14	AUR4514
Ø4.5X16	AUR4516
Ø4.5X18	AUR4518
Ø4.5X20	AUR4520



## AURA BOLT SCREW

Size	Reference Number
H12	AEXPBS12
H14	AEXPBS14
H16	AEXPBS16
H18	AEXPBS18
H20	AEXPBS20



SF-D-3-3/24.08.2020/R00

### Address

Mustafa Kemal Mahallesi, 2133. Sok No: 4/2  
Royale Office PC: 06510 Çankaya/Ankara, Turkey

### Phone

+90 312 473 82 80

### Email

sales@osimplant.com.tr