

JUVE PEDIATRIC SPINAL FIXATION SYSTEM



**Perfect for
pediatric
spine
stabilization**

JUVE Pediatric Spinal Fixation System HIGHLIGHTS

- ♀ Smallest Cup sizes compared to other screws used in the same area
- ♀ Compatible with 4.5mm rod, Jube Polyaxial and monoaxial screws are starting with 20mm and goes up to 75mm with 5mm increments
- ♀ Similar systems are %27 bigger
- ♀ Tav, Tan & CoCr Rod materials options
- ♀ Compatible with 3.00mm Porthos posterior cervical fixation system





JUVE PEDIATRIC SPINAL FIXATION SYSTEM

PURPOSE OF USAGE

JUVE Non- sterile Pediatric Spinal Fixation Systems aim to return the segmental movement and load pattern of the spine to normal, and to reduce the load on the disk and posterior elements. They are spinal fixation systems that can be used with hooks, rods and connectors in degenerative spinal diseases and spinal deformities.

INDICATIONS

- ♀ Degenerative spinal disease
- ♀ Segmental instability due to degenerative disc disease
- ♀ Spondylolisthesis and spondylolysis
- ♀ Trauma (slip, dislocation, etc.)
- ♀ Spinal stenosis (narrow canal)
- ♀ Scoliosis
- ♀ Kyphos
- ♀ Lordosis
- ♀ Scheuermann's hunchback
- ♀ Tumor, infection
- ♀ Pseudoarthrosis
- ♀ Spine fractures
- ♀ Relapse disc hernia





JUVE PEDIATRIC SPINAL FIXATION SYSTEM



JUVE POLYAXIAL SCREW

Size	Reference Number
Ø4,0X20	JPOL4020
Ø4,0X25	JPOL4025
Ø4,0X30	JPOL4030
Ø4,0X35	JPOL4035
Ø4,0X40	JPOL4040
Ø4,5X20	JPOL4520
Ø4,5X25	JPOL4525
Ø4,5X30	JPOL4530
Ø4,5X35	JPOL4535
Ø4,5X40	JPOL4540



Size	Reference Number
Ø5,0X20	JPOL5020
Ø5,0X25	JPOL5025
Ø5,0X30	JPOL5030
Ø5,0X35	JPOL5035
Ø5,0X40	JPOL5040
Ø5,0X45	JPOL5045
Ø5,0X50	JPOL5050
Ø5,5X20	JPOL5520
Ø5,5X25	JPOL5525
Ø5,5X30	JPOL5530
Ø5,5X35	JPOL5535
Ø5,5X40	JPOL5540
Ø5,5X45	JPOL5545
Ø5,5X50	JPOL5550

Size	Reference Number
Ø6,0X30	JPOL6030
Ø6,0X35	JPOL6035
Ø6,0X40	JPOL6040
Ø6,0X45	JPOL6045
Ø6,0X50	JPOL6050
Ø6,0X55	JPOL6055
Ø6,0X60	JPOL6060
Ø6,0X65	JPOL6065
Ø6,0X70	JPOL6070
Ø6,5X75	JPOL6575



JUVE PEDIATRIC SPINAL FIXATION SYSTEM



JUVE MONOAXIAL SCREW

Size	Reference Number
Ø4,0X20	JMON4020
Ø4,0X25	JMON4025
Ø4,0X30	JMON4030
Ø4,0X35	JMON4035
Ø4,0X40	JMON4040
Ø4,5X20	JMON4520
Ø4,5X25	JMON4525
Ø4,5X30	JMON4530
Ø4,5X35	JMON4535
Ø4,5X40	JMON4540

Size	Reference Number
Ø5,0X20	JMON5020
Ø5,0X25	JMON5025
Ø5,0X30	JMON5030
Ø5,0X35	JMON5035
Ø5,0X40	JMON5040
Ø5,0X45	JMON5045
Ø5,0X50	JMON5050
Ø5,5X20	JMON5520
Ø5,5X25	JMON5525
Ø5,5X30	JMON5530
Ø5,5X35	JMON5535
Ø5,5X40	JMON5540
Ø5,5X45	JMON5545
Ø5,5X50	JMON5550




Size	Reference Number
Ø6,0X30	JMON6030
Ø6,0X35	JMON6035
Ø6,0X40	JMON6040
Ø6,0X45	JMON6045
Ø6,0X50	JMON6050
Ø6,0X55	JMON6055
Ø6,0X60	JMON6060
Ø6,0X65	JMON6065
Ø6,0X70	JMON6070
Ø6,0X75	JMON6075

JUVE PEDIATRIC SPINAL FIXATION SYSTEM



JUVE POLYAXIAL REDUCTION SCREW


Size	Reference Number
Ø4,0X20	JPOLRED4020
Ø4,0X25	JPOLRED4025
Ø4,0X30	JPOLRED4030
Ø4,0X35	JPOLRED4035
Ø4,0X40	JPOLRED4040
Ø4,5X20	JPOLRED4520
Ø4,5X25	JPOLRED4525
Ø4,5X30	JPOLRED4530
Ø4,5X35	JPOLRED4535
Ø4,5X40	JPOLRED4540
Ø4,5X45	JPOLRED4545



Size	Reference Number
Ø5,0X20	JPOLRED5020
Ø5,0X25	JPOLRED5025
Ø5,0X30	JPOLRED5030
Ø5,0X35	JPOLRED5035
Ø5,0X40	JPOLRED5040
Ø5,0X45	JPOLRED5045
Ø5,0X50	JPOLRED5050
Ø5,5X20	JPOLRED5520
Ø5,5X25	JPOLRED5525
Ø5,5X30	JPOLRED5530
Ø5,5X35	JPOLRED5535
Ø5,5X40	JPOLRED5540
Ø5,5X45	JPOLRED5545
Ø5,5X50	JPOLRED5550



Size	Reference Number
Ø6,0X30	JPOLRED6030
Ø6,0X35	JPOLRED6035
Ø6,0X40	JPOLRED6040
Ø6,0X45	JPOLRED6045
Ø6,0X50	JPOLRED6050
Ø6,0X55	JPOLRED6055
Ø6,0X60	JPOLRED6060
Ø6,0X65	JPOLRED6065
Ø6,0X70	JPOLRED6070
Ø6,0X75	JPOLRED6075



JUVE PEDIATRIC SPINAL FIXATION SYSTEM



JUVE MONOAXIAL REDUCTION SCREW

Size	Reference Number
Ø4,0X20	JMONRED4020
Ø4,0X25	JMONRED4025
Ø4,0X30	JMONRED4030
Ø4,0X35	JMONRED4035
Ø4,0X40	JMONRED4040
Ø4,5X20	JMONRED4520
Ø4,5X25	JMONRED4525
Ø4,5X30	JMONRED4530
Ø4,5X35	JMONRED4535
Ø4,5X40	JMONRED4540

Size	Reference Number
Ø5,0X20	JMONRED5020
Ø5,0X25	JMONRED5025
Ø5,0X30	JMONRED5030
Ø5,0X35	JMONRED5035
Ø5,0X40	JMONRED5040
Ø5,0X45	JMONRED5045
Ø5,0X50	JMONRED5050



Size	Reference Number
Ø5,5X20	JMONRED5520
Ø5,5X25	JMONRED5525
Ø5,5X30	JMONRED5530
Ø5,5X35	JMONRED5535
Ø5,5X40	JMONRED5540
Ø5,5X45	JMONRED5545
Ø5,5X50	JMONRED5550

Size	Reference Number
Ø6,0X30	JMONRED6030
Ø6,0X35	JMONRED6035
Ø6,0X40	JMONRED6040
Ø6,0X45	JMONRED6045
Ø6,0X50	JMONRED6050
Ø6,0X55	JMONRED6055
Ø6,0X60	JMONRED6060
Ø6,0X65	JMONRED6065
Ø6,0X70	JMONRED6070
Ø6,0X75	JMONRED6075



JUVE PEDIATRIC SPINAL FIXATION SYSTEM

JUVE ROD (Ti6Al4V Eli)-(TAV)

Size	Reference Number
Ø4.5X60	JRTAV60
Ø4.5X80	JRTAV80
Ø4.5X100	JRTAV100
Ø4.5X120	JRTAV120
Ø4.5X150	JRTAV150
Ø4.5X200	JRTAV200
Ø4.5X400	JRTAV400
Ø4.5X500	JRTAN45500
Ø4.5X500	JRTAV45500
Ø4.5X500	JRCOCR45500



JUVE OPEN LATERAL CONNECTOR

Size	Reference Number
10mm	JOLC10
15mm	JOPNC15
20mm	JOLC20
25mm	JOLC25
30mm	JOLC30
40mm	JOPNC40
50mm	JJLC50



JUVE HOOK CONNECTOR

Size	Reference Number
30mm	JHC30
40mm	JHC40
50mm	JHC50
60mm	JHC60
70mm	JHC70
80mm	JHC80



JUVE AXIAL CONNECTOR

Size	Reference Number
65mm	JCAC65
85mm	JCAC85
110mm	JCAC110



JUVE PEDIATRIC SPINAL FIXATION SYSTEM

JUVE CLOSED LATERAL CONNECTOR

Size	Reference Number
10mm	JCLC10
15mm	JCLC15
20mm	JCLC20
25mm	JCLC25
30mm	JOLC30
40mm	JOLC40
50mm	JOLC50



JUVE DOUBLE CUP CONNECTOR

Size	Reference Number
S	JDCDS
M	JDCDM
L	JDCDL



JUVE DOMINO CONNECTOR - SHORT

Size	Reference Number
Ø4.5XØ4.5	JDCMLCS4646
Ø4.5XØ5.5	JDCMLCS4656
Ø4.5XØ6.0	JDCMLCS4661
Ø4.5XØ3.0	JDCMLCS4631



JUVE DOMINO CONNECTOR - SHORT

Size	Reference Number
Ø4.5XØ4.5	JDCMLCL4646
Ø4.5XØ5.5	JDCMLCL4656
Ø4.5XØ6.0	JDCMLCL4661
Ø4.5XØ3.0	JDCMLCL4631





JUVE PEDIATRIC SPINAL FIXATION SYSTEM

JUVE DOMINO CONNECTOR - SHORT STRAIGHT

Size	Reference Number
Ø4.5XØ4.5	JDCSCL4646
Ø4.5XØ5.5	JDCSCL4656
Ø4.5XØ6.0	JDCSCL4661
Ø4.5XØ3.0	JDCSCL4631



JUVE DOMINO CONNECTOR - SHORT STRAIGHT

Size	Reference Number
Ø4.5XØ4.5	JDCSCS4646
Ø4.5XØ5.5	JDCSCS4656
Ø4.5XØ6.0	JDCSCS4661
Ø4.5XØ3.0	JDCSCS4631




JUVE LAMINAR ANGLED BLADE HOOK

Size	Reference Number
STANDARD SIZE	JLABH




JUVE PEDICLE HOOK

Size	Reference Number
STANDARD SIZE	JPEHOS



JUVE EXTENDED HOOK

Size	Reference Number
STANDARD SIZE	JEH





JUVE PEDIATRIC SPINAL FIXATION SYSTEM

JUVE LAMINAR HOOK NARROW BLADE

Size	Reference Number
S	JLHNB-S
M	JLHNB-M
L	JLHNB-L




JUVE LAMINAR HOOK WIDE BLADE

Size	Reference Number
S	JLHWB-S
M	JLHWB-M
L	JLHWB-L




JUVE ANGLED HOOK NARROW BLADE-RIGHT

Size	Reference Number
S	JAHNB-SR
M	JAHNB-MR
L	JAHNB-LR




JUVE ANGLED HOOK NARROW BLADE-LEFT

Size	Reference Number
S	JAHNB-SL
M	JAHNB-ML
L	JAHNB-LL




JUVE OFFSET HOOK NARROW BLADE-RIGHT

Size	Reference Number
S	JOHNB-SR
M	JOHNB-MR
L	JOHNB-LR




JUVE OFFSET HOOK NARROW BLADE-LEFT

Size	Reference Number
S	JOHNB-SL
M	JOHNB-ML
L	JOHNB-LL



JUVE NUT



Size	Reference Number
S	JOHNB-SR

SF-G-3-3/24.08.2020/R00

Address

Mustafa Kemal Mahallesi, 2133. Sok No: 4/2
Royale Office PC: 06510 Çankaya/Ankara, Turkey

Phone

+90 312 473 82 80

Email

sales@osimplant.com.tr