

ON PLUS SPINAL FIXATION SYSTEM



**Perfect for
spine
stabilization**

ON PLUS Spinal Fixation System HIGHLIGHTS

- ♀ Titanium raw material
- ♀ Ideal for use when a larger rod diameter is needed
- ♀ Provides technological advancements to treat a wide range of spinal pathologies
- ♀ Offers a comprehensive solution for rigid posterior fixation of the thoracolumbar regions of the spine
- ♀ Strong structure 6.mm rod
- ♀ Simple and easy usage
- ♀ Wide range of implant choices
- ♀ Screws are designed with DUAL LEAD-DUAL CORE technique with specific threaded body






ON PLUS SPINAL FIXATION SYSTEM

PURPOSE OF USAGE

ON PLUS Non-Sterile Spinal Fixation System is indicated to provide spinal alignment and stabilization. ON PLUS Non-Sterile Spinal Fixation System provides mechanical and biological fixation for early postoperative periods, and provide rigid, rapid and effective fusion. They support the spine during recovery. It can be applied to one or more levels according to various characteristics of the patient and the disease. The components are available in a variety of lengths in order to accommodate patient anatomy. The implants are compatible with all 6 mm rod systems. Screws are designed Dual Lead-Dual Core technique with specific threaded body. Smart screw design eliminates need to tap.

INDICATIONS

- ♀ Degenerative spinal diseases
- ♀ Segmental instability due to degenerative disc disease
- ♀ Spondylolisthesis and spondylolysis
- ♀ Trauma (slip, dislocation, etc.)
- ♀ Spinal stenosis (narrow canal)
- ♀ Scoliosis
- ♀ Kyphosis
- ♀ Lordosis
- ♀ Scheuermann's hunchback
- ♀ Tumor, infection





ON PLUS SPINAL FIXATION SYSTEM

ON PLUS POLYAXIAL SCREW

Size	Reference Number
Ø4,5X20	ON0102-4520
Ø4,5X25	ON0102-4525
Ø4,5X30	ON0102-4530
Ø4,5X35	ON0102-4535
Ø4,5X40	ON0102-4540
Ø4,5X45	ON0102-4545
Ø4,5X50	ON0102-4550
Ø4,5X55	ON0102-4555
Ø4,5X60	ON0102-4560
Ø4,5X65	ON0102-4565
Ø4,5X70	ON0102-4570
Ø4,5X75	ON0102-4575
Ø4,5X80	ON0102-4580
Ø4,5X85	ON0102-4585
Ø4,5X90	ON0102-4590
Ø4,5X95	ON0102-4595
Ø4,5X100	ON0102-45100



Size	Reference Number
Ø6,5X20	ON0102-6520
Ø6,5X25	ON0102-6525
Ø6,5X30	ON0102-6530
Ø6,5X35	ON0102-6535
Ø6,5X40	ON0102-6540
Ø6,5X45	ON0102-6545
Ø6,5X50	ON0102-6550
Ø6,5X55	ON0102-6555
Ø6,5X60	ON0102-6560
Ø6,5X65	ON0102-6565
Ø6,5X70	ON0102-6570
Ø6,5X75	ON0102-6575
Ø6,5X80	ON0102-6580
Ø6,5X85	ON0102-6585
Ø6,5X90	ON0102-6590
Ø6,5X95	ON0102-6595
Ø6,5X100	ON0102-65100



Size	Reference Number
Ø5,5X20	ON0102-5520
Ø5,5X25	ON0102-5525
Ø5,5X30	ON0102-5530
Ø5,5X35	ON0102-5535
Ø5,5X40	ON0102-5540
Ø5,5X45	ON0102-5545
Ø5,5X50	ON0102-5550
Ø5,5X55	ON0102-5555
Ø5,5X60	ON0102-5560
Ø5,5X65	ON0102-5565
Ø5,5X70	ON0102-5570
Ø5,5X75	ON0102-5570
Ø5,5X80	ON0102-5580
Ø5,5X85	ON0102-5585
Ø5,5X90	ON0102-5590
Ø5,5X95	ON0102-5595
Ø5,5X100	ON0102-55100





ON PLUS SPINAL FIXATION SYSTEM

ON PLUS POLYAXIAL SCREW

Size	Reference Number
Ø7,5X20	ON-PR7520
Ø7,5X25	ON-PR7525
Ø7,5X30	ON-PR7530
Ø7,5X35	ON-PR7535
Ø7,5X40	ON-PR7540
Ø7,5X45	ON-PR7545
Ø7,5X50	ON-PR7550
Ø7,5X55	ON-PR7555
Ø7,5X60	ON-PR7560
Ø7,5X65	ON-PR7565
Ø7,5X70	ON-PR7570
Ø7,5X75	ON-PR7575
Ø7,5X80	ON-PR7580
Ø7,5X85	ON-PR7585
Ø7,5X90	ON-PR7590
Ø7,5X95	ON-PR7595
Ø7,5X100	ON-PR75100



Size	Reference Number
Ø8,5X20	ON-PR8520
Ø8,5X25	ON-PR8525
Ø8,5X30	ON-PR8530
Ø8,5X35	ON-PR8535
Ø8,5X40	ON-PR8540
Ø8,5X45	ON-PR8545
Ø8,5X50	ON-PR8550
Ø8,5X55	ON-PR8555
Ø8,5X60	ON-PR8560
Ø8,5X65	ON-PR8565
Ø8,5X70	ON-PR8570
Ø8,5X75	ON-PR8575
Ø8,5X80	ON-PR8580
Ø8,5X85	ON-PR8585
Ø8,5X90	ON-PR8590
Ø8,5X95	ON-PR8595
Ø8,5X100	ON-PR85100

ON PLUS POLYAXIAL REDUCTION SCREW

Size	Reference Number
Ø4,5X20	ON0122-4520
Ø4,5X25	ON0122-4525
Ø4,5X30	ON0122-4530
Ø4,5X35	ON0122-4535
Ø4,5X40	ON0122-4540
Ø4,5X45	ON0122-4545
Ø4,5X50	ON0122-4550
Ø4,5X55	ON0122-4555
Ø4,5X60	ON0122-4560
Ø4,5X65	ON0122-4565
Ø4,5X70	ON0122-4570
Ø4,5X75	ON0122-4575
Ø4,5X80	ON0122-4580
Ø4,5X85	ON0122-4585
Ø4,5X90	ON0122-4590
Ø4,5X95	ON0122-4595
Ø4,5X100	ON0122-45100





ON PLUS POLYAXIAL REDUCTION SCREW

Size	Reference Number
Ø6,5X20	ON0122-6520
Ø6,5X25	ON0122-6525
Ø6,5X30	ON0122-6530
Ø6,5X35	ON0122-6535
Ø6,5X40	ON0122-6540
Ø6,5X45	ON0122-6545
Ø6,5X50	ON0122-6550
Ø6,5X55	ON0122-6555
Ø6,5X60	ON0122-6560
Ø6,5X65	ON0122-6565
Ø6,5X70	ON0122-6570
Ø6,5X75	ON0122-6575
Ø6,5X80	ON0122-6580
Ø6,5X85	ON0122-6585
Ø6,5X90	ON0122-6590
Ø6,5X95	ON0122-6595
Ø6,5X100	ON0122-65100



Size	Reference Number
Ø7,5X20	ON0122-7520
Ø7,5X25	ON0122-7525
Ø7,5X30	ON0122-7530
Ø7,5X35	ON0122-7535
Ø7,5X40	ON0122-7540
Ø7,5X45	ON0122-7545
Ø7,5X50	ON0122-7550
Ø7,5X55	ON0122-7555
Ø7,5X60	ON0122-7560
Ø7,5X65	ON0122-7565
Ø7,5X70	ON0122-7570
Ø7,5X75	ON0122-7575
Ø7,5X80	ON0122-7580
Ø7,5X85	ON0122-7585
Ø7,5X90	ON0122-7590
Ø7,5X95	ON0122-7595
Ø7,5X100	ON0122-75100

Size	Reference Number
Ø5,5X20	ON0122-5520
Ø5,5X25	ON0122-5525
Ø5,5X30	ON0122-5530
Ø5,5X35	ON0122-5535
Ø5,5X40	ON0122-5540
Ø5,5X45	ON0122-5545
Ø5,5X50	ON0122-5550
Ø5,5X55	ON0122-5555
Ø5,5X60	ON0122-5560
Ø5,5X65	ON0122-5565
Ø5,5X70	ON0122-5570
Ø5,5X75	ON0122-5575
Ø5,5X80	ON0122-5580
Ø5,5X85	ON0122-5585
Ø5,5X90	ON0122-5590
Ø5,5X95	ON0122-5595
Ø5,5X100	ON0122-55100



ON PLUS MONOAXIAL SCREW

Size	Reference Number
Ø4,5X20	ON0101-4520
Ø4,5X25	ON0101-4525
Ø4,5X30	ON0101-4530
Ø4,5X35	ON0101-4535
Ø4,5X40	ON0101-4540
Ø4,5X45	ON0101-4545
Ø4,5X50	ON0101-4550
Ø4,5X55	ON0101-4555
Ø4,5X60	ON0101-4560
Ø5,5X20	ON0101-5520
Ø5,5X25	ON0101-5525
Ø5,5X30	ON0101-5530
Ø5,5X35	ON0101-5535
Ø5,5X40	ON0101-5540
Ø5,5X45	ON0101-5545
Ø5,5X50	ON0101-5550
Ø5,5X55	ON0101-5555
Ø5,5X60	ON0101-5560
Ø6,5X20	ON0101-6520
Ø6,5X25	ON0101-6525
Ø6,5X30	ON0101-6530
Ø6,5X35	ON0101-6535
Ø6,5X40	ON0101-6540
Ø6,5X45	ON0101-6545
Ø6,5X50	ON0101-6550
Ø6,5X55	ON0101-6555
Ø6,5X60	ON0101-6560

Size	Reference Number
Ø7,5X20	ON0101-7520
Ø7,5X25	ON0101-7525
Ø7,5X30	ON0101-7530
Ø7,5X35	ON0101-7535
Ø7,5X40	ON0101-7540
Ø7,5X45	ON0101-7545
Ø7,5X50	ON0101-7550
Ø7,5X55	ON0101-7555
Ø7,5X60	ON0101-7560
Ø8,5X20	ON0101-8520
Ø8,5X25	ON0101-8525
Ø8,5X30	ON0101-8530
Ø8,5X35	ON0101-8535
Ø8,5X40	ON0101-8540
Ø8,5X45	ON0101-8545
Ø8,5X50	ON0101-8550
Ø8,5X55	ON0101-8555
Ø8,5X60	ON0101-8560





ON PLUS MONOAXIAL REDUCTION SCREW

Size	Reference Number
Ø4,5X20	ON0103-4520
Ø4,5X25	ON0103-4525
Ø4,5X30	ON0103-4530
Ø4,5X35	ON0103-4535
Ø4,5X45	ON0103-4545
Ø4,5X50	ON0103-4550
Ø4,5X55	ON0103-4555
Ø4,5X60	ON0103-4560
Ø5,5X20	ON0103-5520
Ø5,5X25	ON0103-5525
Ø5,5X30	ON0103-5530
Ø5,5X35	ON0103-5535
Ø5,5X40	ON0103-5540
Ø5,5X45	ON0103-5545
Ø5,5X50	ON0103-5550
Ø5,5X55	ON0103-5555
Ø5,5X60	ON0103-5560
Ø6,5X20	ON0103-6520
Ø6,5X25	ON0103-6525
Ø6,5X30	ON0103-6530
Ø6,5X35	ON0103-6535
Ø6,5X40	ON0103-6540
Ø6,5X45	ON0103-6545
Ø6,5X50	ON0103-6550
Ø6,5X55	ON0103-6555
Ø6,5X60	ON0103-6560

Size	Reference Number
Ø7,5X20	ON0103-7520
Ø7,5X25	ON0103-7525
Ø7,5X30	ON0103-7530
Ø7,5X35	ON0103-7535
Ø7,5X40	ON0103-7540
Ø7,5X45	ON0103-7545
Ø7,5X50	ON0103-7550
Ø7,5X55	ON0103-7555
Ø7,5X60	ON0103-7560
Ø8,5X20	ON0103-8520
Ø8,5X25	ON0103-8525
Ø8,5X30	ON0103-8530
Ø8,5X35	ON0103-8535
Ø8,5X40	ON0103-8540
Ø8,5X45	ON0103-8545
Ø8,5X50	ON0103-8550
Ø8,5X55	ON0103-8555
Ø8,5X60	ON0103-8560



ON PLUS SPINAL FIXATION SYSTEM



ON PLUS ROD

Size	Reference Number
Ø6.0X40	ON0150-0040
Ø6.0X45	ON0150-0045
Ø6.0X50	ON0150-0050
Ø6.0X60	ON0150-0060
Ø6.0X70	ON0150-0070
Ø6.0X80	ON0150-0080
Ø6.0X90	ON0150-0090
Ø6.0X100	ON0150-0100
Ø6.0X110	ON0150-0110
Ø6.0X120	ON0150-0120
Ø6.0X130	ON0150-0130
Ø6.0X140	ON0150-0140
Ø6.0X150	ON0150-0150
Ø6.0X160	ON0150-0160
Ø6.0X170	ON0150-0170
Ø6.0X180	ON0150-0180
Ø6.0X190	ON0150-0190
Ø6.0X200	ON0150-0200
Ø6.0X210	ON0150-0210
Ø6.0X220	ON0150-0220
Ø6.0X230	ON0150-0230
Ø6.0X240	ON0150-0240
Ø6.0X250	ON0150-0250
Ø6.0X260	ON0150-0260
Ø6.0X270	ON0150-0270
Ø6.0X280	ON0150-0280
Ø6.0X290	ON0150-0290
Ø6.0X300	ON0150-0300
Ø6.0X350	ON0150-0350
Ø6.0X400	ON0150-0400
Ø6.0X450	ON0150-0450
Ø6.0X500	ON0150-0500
Ø6.0X600	ON0150-0600
Ø6.0X700	ON0150-0700

Ø6.0x100mm



ON PLUS ROD (Co.Cr)

Size	Reference Number
Ø6.0X40	COCR-040
Ø6.0X50	COCR-050
Ø6.0X60	COCR-060
Ø6.0X70	COCR-070
Ø6.0X80	COCR-080
Ø6.0X90	COCR-090
Ø6.0X100	COCR-100
Ø6.0X110	COCR-110
Ø6.0X120	COCR-120
Ø6.0X130	COCR-130
Ø6.0X140	COCR-140
Ø6.0X150	COCR-150
Ø6.0X160	COCR-160
Ø6.0X170	COCR-170
Ø6.0X180	COCR-180
Ø6.0X190	COCR-190
Ø6.0X200	COCR-200
Ø6.0X210	COCR-210
Ø6.0X220	COCR-220
Ø6.0X230	COCR-230
Ø6.0X240	COCR-240
Ø6.0X250	COCR-250
Ø6.0X260	COCR-260
Ø6.0X270	COCR-270
Ø6.0X280	COCR-280
Ø6.0X290	COCR-290
Ø6.0X300	COCR-300
Ø6.0X350	COCR-350
Ø6.0X400	COCR-400
Ø6.0X450	COCR-450
Ø6.0X500	COCR-500
Ø6.0X600	COCR-600
Ø6.0X700	COCR-700

Ø6.0x100mm



ON PLUS SPINAL FIXATION SYSTEM



ON PLUS CURVED ROD

Size	Reference Number
Ø6.0X40	ONCR040
Ø6.0X50	ONCR050
Ø6.0X60	ONCR060
Ø6.0X70	ONCR070
Ø6.0X80	ONCR080
Ø6.0X90	ONCR090
Ø6.0X100	ONCR0100
Ø6.0X110	ONCR0110
Ø6.0X120	ONCR0120
Ø6.0X130	ONCR0130
Ø6.0X140	ONCR0140
Ø6.0X150	ONCR0150
Ø6.0X160	ONCR0160
Ø6.0X170	ONCR0170
Ø6.0X180	ONCR0180
Ø6.0X190	ONCR0190
Ø6.0X200	ONCR0200



ON PLUS DYNAMIC PEEK ROD

Size	Reference Number
Ø6.0X30	ONPDPR30
Ø6.0X40	ONPDPR40
Ø6.0X50	ONPDPR50
Ø6.0X60	ONPDPR60
Ø6.0X70	ONPDPR70
Ø6.0X80	ONPDPR80
Ø6.0X90	ONPDPR90
Ø6.0X100	ONPDPR100
Ø6.0X110	ONPDPR110
Ø6.0X120	ONPDPR120
Ø6.0X130	ONPDPR130
Ø6.0X135	ONPDPR135
Ø6.0X140	ONPDPR140



ON PLUS STRAIGHT PEEK ROD

Size	Reference Number
Ø6.0X100	ONSPR100
Ø6.0X200	ONSPR200

Ø6.0x100mm



ON PLUS DYNAMIC ROD

Size	Reference Number
Ø6.0X60	HER-60
Ø6.0X70	HER-70
Ø6.0X80	HER-80
Ø6.0X90	HER-90
Ø6.0X100	HER-100
Ø6.0X110	HER-110
Ø6.0X120	HER-120
Ø6.0X130	HER-130
Ø6.0X140	HER-140





ON PLUS SPINAL FIXATION SYSTEM


ON PLUS FLAT ROD

Size	Reference Number
Ø6X40	F0040
Ø6X50	F0050
Ø6X60	F0060
Ø6X70	F0070
Ø6X80	F0080
Ø6X90	F0090
Ø6X100	F0100
Ø6X110	F0110
Ø6X120	F0120
Ø6X130	F0130
Ø6X140	F0140
Ø6X150	F0150
Ø6X160	F0160
Ø6X170	F0170
Ø6X180	F0180
Ø6X190	F0190
Ø6X200	F0200
Ø6X250	F0250
Ø6X300	F0300
Ø6X400	F0400
Ø6X500	F0500
Ø6X600	F0600
Ø6X700	F0700



ON PLUS MULTIAXIAL TRANSVERSE CONNECTOR

Size	Reference Number
33-40	ON0120-0035
40-50	ON0120-0038
50-60	ON0120-0050
60-70	ON0120-0042



ON PLUS HOOK TRANSVERSE CONNECTOR

Size	Reference Number
50mm	ONHTC50
60mm	ONHTC60
70mm	ONHTC70
80mm	ONHTC80
90mm	ONHTC90
100mm	ONHTC100






ON PLUS SPINAL FIXATION SYSTEM


ON PLUS LATERAL CONNECTOR(CLOSED)

Size	Reference Number
20mm	ON0192-0620
25mm	ON0192-0625
30mm	ON-LCC30
35mm	ON-LCC35
40mm	ON-LCC40
45mm	ON-LCC45
50mm	ON-LCC50
55mm	ON-LCC55
60mm	ON-LCC60
65mm	ON-LCC65
70mm	ON-LCC70



ON PLUS LATERAL CONNECTOR(OPENED)

Size	Reference Number
20mm	ON-LC020
25mm	ON-LC025
30mm	ON-LC030
35mm	ON-LC035
40mm	ON-LC040



ON PLUS DOUBLE CUP DOMINO



Size	Reference Number
S	DOM-ONDSCS
M	DOM-ONDSCM
L	DOM-ONDSCCL

ON PLUS NUT



Size	Reference Number
T30	ONPLUSNUT

ON PLUS SPINAL FIXATION SYSTEM



ON PLUS LAMINA HOOK - NARROW BLADE

Size	Reference Number
SMALL	ONPLHNB-S
MEDIUM	ONPLHNB-M
LARGE	ONPLHNB-L
XLARGE	ONPLHNB-XL



ON PLUS LAMINA HOOK - WIDE BLADE

Size	Reference Number
SMALL	ONPLHNB-S
MEDIUM	ONPLHNB-M
LARGE	ONPLHNB-L
XLARGE	ONPLHNB-XL

ON PLUS OFFSET LAMINA HOOK NARROW BLADE - LEFT

Size	Reference Number
SMALL	ONPOHNB-SL
MEDIUM	ONPOHNB-ML
LARGE	ONPOHNB-LL
XLARGE	ONP-LHOL0511

ON PLUS OFFSET LAMINA HOOK NARROW BLADE - RIGHT

Size	Reference Number
SMALL	ONPOHNB-SR
MEDIUM	ONPOHNB-MR
LARGE	ONPOHNB-LR
XLARGE	ONP-LHOR0511

ON PLUS ANGLED HOOK NARROW BLADE LEFT


Size	Reference Number
SMALL	ONP-LHAL0505
MEDIUM	ONP-LHAL0507
LARGE	ONP-LHAL0509

ON PLUS ANGLED HOOK NARROW BLADE RIGHT

Size	Reference Number
SMALL	ONP-LHAR0505
MEDIUM	ONP-LHAR0507
LARGE	ONP-LHAR0509




ON PLUS PEDICLE HOOK




Size	Reference Number
STANDARD	ONP-PH0607

ON PLUS EXTENDED HOOK



Size	Reference Number
STANDARD	ONPSEH

ON PLUS LAMINAR ANGLED BLADE HOOK



Size	Reference Number
STANDARD	ONPLABH