



# TITANOPEEK ACF STAND ALONE CERVICAL CAGE SYSTEM




**Perfect for  
provide  
cervical  
fusion.**

## **TITANOPEEK ACF Stand Alone Cervical Cage System HIGHLIGHTS**

- ☞ Integrate a hollow PEEK cage with a titanium screw locking mechanism
- ☞ Can be safely applied to more than one level
- ☞ Large bone graft area for better fusion
- ☞ No need for anterior profile
- ☞ Self-drilling screws , Self-locking screws, Stand alone cage
- ☞ Offer in a variety of sizes



# TITANOPEEK ACF STAND ALONE CERVICAL CAGE SYSTEM




## PURPOSE OF USAGE

TITANOPEEK ACF STAND ALONE CERVICAL CAGE SYSTEM, which are made of biocompatible polyetheretherketone (PEEK) raw material, the rest of which are made of Ti6Al4V alloy raw material, restores the disc distance, increases the fusion area, provides strong end-plate adhesion thanks to the threaded structure and increases the fusion area. It supports bone fusion and shortens the operative time. Tantalum-made marker pins allow the cage to be displayed on the X-Ray.


Our TITANOPEEK cage products are used with screws. Screws are produced from Ti6Al4V alloy raw material

## INDICATIONS

- ♀ Vertebral body destruction with kyphosis
- ♀ Segmental instability
- ♀ Epidural / paravertebral abscess requiring severe bone resection
- ♀ Stabilization after treatment of a herniated cervical intervertebral disc
- ♀ Osteophytosis compressing the nerve roots and/or spinal cord
- ♀ Degenerative intervertebral instability
- ♀ Degenerative disc diseases and spinal instability
- ♀ Cervical spine trauma
- ♀ Radiculopathy
- ♀ Myelopathy
- ♀ Osteomyelitis




# TITANOPEEK ACF STAND ALONE CERVICAL CAGE SYSTEM



## TITANOPEEK HYBRID CAGE


Size	Reference Number
12X14 H5 - 0°	TNO-1214H5
12X14 H6 - 0°	TNO-1214H6
12X14 H7 - 0°	TNO-1214H7
12X14 H8 - 0°	TNO-1214H8
12X14 H5 - 5° LORDOSIS	TNL-1214H5
12X14 H6 - 5° LORDOSIS	TNL-1214H6
12X14 H7 - 5° LORDOSIS	TNL-1214H7
12X14 H8 - 5° LORDOSIS	TNL-1214H8

Size	Reference Number
13X15 H5 - 0°	TNO-1315H5
13X15 H6 - 0°	TNO-1315H6
13X15 H7 - 0°	TNO-1315H7
13X15 H8 - 0°	TNO-1315H8
13X15 H5 - 5° LORDOSIS	TNL-1315H5
13X15 H6 - 5° LORDOSIS	TNL-1315H6
13X15 H7 - 5° LORDOSIS	TNL-1315H7
13X15 H8 - 5° LORDOSIS	TNL-1315H8



Size	Reference Number
14X16 H5 - 0°	TNO-1416H5
14X16 H6 - 0°	TNO-1416H6
14X16 H7 - 0°	TNO-1416H7
14X16 H8 - 0°	TNO-1416H8
14X16 H5 - 5° LORDOSIS	TNL-1416H5
14X16 H6 - 5° LORDOSIS	TNL-1416H6
14X16 H7 - 5° LORDOSIS	TNL-1416H7
14X16 H8 - 5° LORDOSIS	TNL-1416H8



# TITANOPEEK ACF STAND ALONE CERVICAL CAGE SYSTEM



## TITANOPEEK TITANIUM CAGE


Size	Reference Number
12X14 H5 - 0°	TTNO-1214H5
12X14 H6 - 0°	TTNO-1214H6
12X14 H7 - 0°	TTNO-1214H7
12X14 H8 - 0°	TTNO-1214H8
12X14 H5 - 5° LORDOSIS	TTNL-1214H5
12X14 H6 - 5° LORDOSIS	TTNL-1214H6
12X14 H7 - 5° LORDOSIS	TTNL-1214H7
12X14 H8 - 5° LORDOSIS	TTNL-1214H8

Size	Reference Number
13X15 H5 - 0°	TTNO-1315H5
13X15 H6 - 0°	TTNO-1315H6
13X15 H7 - 0°	TTNO-1315H7
13X15 H8 - 0°	TTNO-1315H8
13X15 H5 - 5° LORDOSIS	TTNL-1315H5
13X15 H6 - 5° LORDOSIS	TTNL-1315H6
13X15 H7 - 5° LORDOSIS	TTNL-1315H7
13X15 H8 - 5° LORDOSIS	TTNL-1315H8



Size	Reference Number
14X16 H5 - 0°	TTNO-1416H5
14X16 H6 - 0°	TTNO-1416H6
14X16 H7 - 0°	TTNO-1416H7
14X16 H8 - 0°	TTNO-1416H8
14X16 H5 - 5° LORDOSIS	TTNL-1416H5
14X16 H6 - 5° LORDOSIS	TTNL-1416H6
14X16 H7 - 5° LORDOSIS	TTNL-1416H7
14X16 H8 - 5° LORDOSIS	TTNL-1416H8

# TITANOPEEK ACF STAND ALONE CERVICAL CAGE SYSTEM



## TITANOPEEK LOCKING SCREW SELF DRILLING

Size	Reference Number
Ø3.5X12	TNSD-3512
Ø3.5X14	TNSD-3514
Ø3.5X16	TNSD-3516
Ø3.5X18	TNSD-3518
Ø3.5X20	TNSD-3520
Ø4.0X12	TNSD-4012
Ø4.0X14	TNSD-4014
Ø4.0X16	TNSD-4016
Ø4.0X18	TNSD-4018
Ø4.0X20	TNSD-4020



## TITANOPEEK LOCKING SCREW BLUNT TYPE



Size	Reference Number
Ø3.5X12	TNBT-3512
Ø3.5X14	TNBT-3514
Ø3.5X16	TNBT-3516
Ø3.5X18	TNBT-3518
Ø3.5X20	TNBT-3520
Ø4.0X12	TNBT-4012
Ø4.0X14	TNBT-4014
Ø4.0X16	TNBT-4016
Ø4.0X18	TNBT-4018
Ø4.0X20	TNBT-4020

SK-A-4-3-3/24.08.2020/R00

**Address**

Mustafa Kemal Mahallesi, 2133. Sok No: 4/2  
Royale Office PC: 06510 Çankaya/Ankara, Turkey

**Phone**

+90 312 473 82 80

**Email**

sales@osimplant.com.tr