


ZELOS PEEK/TITANIUM PLIF CAGE




**Perfect for
provide
interbody
fusion.**

ZELOS PEEK/TITANIUM PLIF CAGE HIGHLIGHTS

- ⊕ Titanium or PEEK raw material
- ⊕ Tantalum X-ray markers
- ⊕ Large segmental bone graft area for better fusion
- ⊕ Ability to treat collapsed disc spaces
- ⊕ Allows extraordinary grip due to threaded structure
- ⊕ Simple and easy usage
- ⊕ Anatomical design
- ⊕ Special machined toothed surfaces for good adhesion






PURPOSE OF USAGE


♀ ZELOS PLIF (Posterior Lumbar Interbody Fusion) Spinal Cages are used in the treatment of low back and lower extremity pain, numbness and partial mobility loss caused by nerve compression caused by disc degeneration between the vertebrae. In this treatment the intervertebral disc is removed, and grafts and / or intervertebral fusion cages are placed in place to fuse adjacent vertebrae. With this treatment, the degenerated intervertebral disc (possible pain source) is eliminated and the height of the disc and anatomical lordosis angle are restored.

INDICATIONS

- ♀ Vertebral body destruction accompanied by kyphosis
- ♀ Segmental instability
- ♀ Epidural / paravertebral abscesses requiring severe bone resection
- ♀ Stabilization after herniated lumbar intervertebral disc treatment
- ♀ Osteophytosis that compresses nerve roots and / or spinal cord
- ♀ Degenerative intervertebral instability
- ♀ Degenerative disc diseases and spinal instabilities
- ♀ Lumbar spine trauma
- ♀ Radiculopathy
- ♀ Myelopathy




ZELOS PEEK/TITANIUM PLIF CAGE




ZELOS PEEK PLIF CAGE (0°)

Size	Reference Number
24X9-H7 - 0°	ZE240907
24X9-H8 - 0°	ZE240908
24X9-H9 - 0°	ZE240909
24X9-H10 - 0°	ZE240910
24X9-H11 - 0°	ZE240911
24X9-H12 - 0°	ZE240912
24X9-H13 - 0°	ZE240913
24X9-H14 - 0°	ZE240914
24X9-H15 - 0°	ZE240915
24X9-H16 - 0°	ZE240916





ZELOS PEEK PLIF CAGE (4° LORDOSIS)

Size	Reference Number
24X9-H7 - 4°	ZE240407
24X9-H8 - 4°	ZE240408
24X9-H9 - 4°	ZE240409
24X9-H10 - 4°	ZE240410
24X9-H11 - 4°	ZE240411
24X9-H12 - 4°	ZE240412
24X9-H13 - 4°	ZE240413
24X9-H14 - 4°	ZE240414
24X9-H15 - 4°	ZE240415
24X9-H16 - 4°	ZE240416



ZELOS PEEK/TITANIUM PLIF CAGE




ZELOS PEEK PLIF CAGE (8° LORDOSIS)

Size	Reference Number
24X9-H7 - 8°	ZE240807
24X9-H8 - 8°	ZE240808
24X9-H9 - 8°	ZE240809
24X9-H10 - 8°	ZE240810
24X9-H11 - 8°	ZE240811
24X9-H12 - 8°	ZE240812
24X9-H13 - 8°	ZE240813
24X9-H14 - 8°	ZE240814
24X9-H15 - 8°	ZE240815
24X9-H16 - 8°	ZE240816

ZELOS TITANIUM PLIF CAGE

Size	Reference Number
24X9-H7 - 0°	TIZLS2409070S
24X9-H8 - 0°	TIZLS2409080S
24X9-H9 - 0°	TIZLS2409090S
24X9-H10 - 0°	TIZLS2409100S
24X9-H11 - 0°	TIZLS2409110S
24X9-H12 - 0°	TIZLS2409120S
24X9-H13 - 0°	TIZLS2409130S



SK-B-4-3-3/24.08.2020/R00

Address

Mustafa Kemal Mahallesi, 2133. Sok No: 4/2
Royale Office PC: 06510 Çankaya/Ankara, Turkey

Phone

+90 312 473 82 80

Email

sales@osimplant.com.tr