


MIS

OPEN



MODULAR Shank and Housing options


VERSATILE MIS to OPEN compatibility

EFFICIENT implant and instrument traceability

317-436-7801
Info@NexxtSpine.com
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
70-095, Rev A


MINIMALLY INVASIVE

MINIMALLY INVASIVE SURGICAL (MIS) PROCEDURES



INTERNAL REDUCTION

MIS Towers feature 20mm of internal reduction.



INNOVATIVE INSTRUMENTS

Rod Inserters lift straight up upon Rod release.



MIS TO OPEN INSTRUMENT COMPATABILITY

TRADITIONAL OPEN

DEGENERATIVE OPEN SURGICAL PROCEDURES



HOUSING OPTIONS

Standard and Reduction Housings available to meet procedure needs.



OPERATIVE FLEXIBILITY

Multiple approaches for in-situ and back table assembly of modular implants.




SHANKS AND HOUSINGS

Shanks available in Ø4.5-11.5mm x 20-140mm long




HOUSING MODULARITY

All Shanks are compatible with all Housings.




IMPLANT TRACEABILITY

Implants feature removable tracking tabs marked with lot information and UDI.



UNIVERSAL HOUSING

Ability to utilize both OPEN and MIS instrumentation for primary and revision procedures.



DUAL DIAMETER SLOT


Allows use of both Ø5.5mm and Ø6.0 Rods.



SHANK OPTIONS

Cannulated / Solid
Uniplanar / Polyaxial
Dual Lead / Corticocancellous

SET SCREWS & RODS



SNAP LOCK

"Click on" attachment feature quickly secures Set Screw to prevent dropping in-situ.



STANDARDIZED T25 DRIVE

All Shanks and Set Screws use the same T25 drive feature.



AUTO RELEASE

Set Screw releases from Set Screw Driver when reaching final locked position.




ROD VERSATILITY

Rods available in straight, curved, titanium, and CoCr configurations.

HIGHLIGHTED INSTRUMENTS

SHANK INSERTERS

Ability to insert Shank and Housing separately or as an assembly.



ROD REDUCER

25mm+ of additional reduction in both MIS and OPEN.



MIS ROD INSERTER

Lordotic or kyphotic placement with "straight up" removal of instrument.




HOUSING ASSEMBLY TOOL

One step procedure to guarantee shank-housing assembly security.




MIS COMPRESSOR DISTRACTOR

Provides both compression and distraction in one powerful instrument.




COMPATIBLE MATRIXX PRODUCTS

MATRIXX® TLIF Oblique




MATRIXX® SA ALIF




MATRIXX® TLIF / rTLIF



MATRIXX® Lateral



Nexxt Spine, LLC
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
(317) 436-7801
Info@NexxtSpine.com
NexxtSpine.com



70-095, Rev A

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit NexxtSpine.com for additional product information.

All rights reserved. All content herein is protected by copyright, trademarks and other intellectual property rights owned by Nexxt Spine, LLC and must not be redistributed, duplicated or disclosed, in whole or in part without the expressed written consent of Nexxt Spine, LLC. This material is intended for health care professionals, the Nexxt Spine sales force and authorized representatives. Distribution to any other recipient is prohibited.