


# Mercury<sup>®</sup> MIS

## Pedicle Screw System



Proven Concept.  
Uniquely Elegant.  
Minimally Invasive.

# Mercury<sup>®</sup> MIS | Pedicle Screw System

The Mercury Classic MIS system is an innovative, simple, and versatile pedicle screw system designed to limit muscle and tissue disruption through a minimally invasive approach while addressing the clinical needs of posterior fixation.

## System Overview

- Tower based delivery system
- Simple loading and unloading of the towers
- Easy to use instrumentation
- One tray add-on to Mercury Classic

## Towers

- Simple loading and unloading
- Large tower slots makes rod insertion easy
- Locking collar ensures secure screw attachment

## Cannulated Screws

- Low profile
- Self-tapping
- 50° total angulation
- Diameters from 5.5mm to 8.5mm
- Lengths from 30mm to 60mm




## Rods

- Bulleted for easier rod insertion
- Lordosed
- 5.5mm rod diameter
- Lengths from 30mm to 120mm



## Reducer

- Easily engages tower and threaded cap inserter for powerful and precise reduction
- 33mm of total reduction capability



## Rod Length Caliper

- Simple and easy to use
- Eliminates the need to guess rod length



# Mercury<sup>®</sup> MIS

## Pedicle Screw System

Screws		
	Diameter	Length
	4.5	30-50
	5.5	30-55
	6.5	30-60
	7.5	30-60
	8.5	30-50

Dimensions expressed in millimeters  
5mm increments.

Lordosed Rods	
Diameter	Length
5.5	30-120

Dimensions expressed in millimeters  
5mm increments up to 60mm, 10mm increments up to 120mm

