


## THORACOLUMBAR PEDICLE SCREW SYSTEM





Phone: (484) 427-7060 cambermedtech.com info@cambermedtech.com




## cambermedtech.com



## THORACOLUMBAR PEDICLE SCREW SYSTEM

- Cannulated screw diameters available
  - 5.5mm, 6.5mm, 7.5mm, 8.5mm\*
- Poly drag in screw head for easy reception of rod
- Blunt-tipped for bi-cortical purchase
- Pre-lordosed\* and Straight rods with self-dissecting bulleted tip
- Easy break screw tabs that leave a low-profile and clean break
- 8.5 N-m (~6.3lbf-ft) Torque locking set screws
  - Quick insertion with self-centering instruments
  - Thread design eliminates cross-threading
- 25mm Threaded Reduction in screw tab
- Controlled, adjustable compression and distraction technique utilizing a fulcrum
- Low-profile Rod Inserter maintains rigid rod control but keeps incision sizes shorter for multi-level cases
- Multiple dilation options available
  - \*Available upon special request with additional lead time







CM-801-001, REV C