


SEURAT™

POSTERIOR LUMBAR
MODULAR MIS FIXATION PLATFORM

SEURAT™ *flex* modular percutaneous polyaxial screws are designed with versatility in mind. Made for true percutaneous implantation, they allow for complete pedicle stabilization, and are designed with a proprietary mid-tower accordion hinge. This system is to be used in conjunction with the PICASSO II™ system.



Flexible tower provides superior usability in crowded areas

Stability ring prevents splay during implantation, reduction and set screw insertion

Laser scored and tapered housing provides a clean break every time

6.0 percutaneous **rod available in straight and curved**, and designed with a bulleted tip, which allows for smooth insertion

25mm of internal threads aid in rod reduction

Dual lead polyaxial screw with an angulation of $\pm 20^\circ$



CTL AMERICA
rethink what's possible®

SEURAT™ Screwbody: 5.5mm Diameter

Part Number	Color	Length (mm)
126.1636	Green	35
126.1641	Green	40
126.1646	Green	45
126.1651	Green	50
126.1656	Green	55

SEURAT™ Flexible MIS Screw Tower

Part Number	Color
126.0034	Blue

SEURAT™ Screwbody: 6.5mm Diameter

Part Number	Color	Length (mm)
126.1836	Gold	35
126.1841	Gold	40
126.1846	Gold	45
126.1851	Gold	50
126.1856	Gold	55
126.1861	Gold	60

SEURAT™ Screwbody: 7.5mm Diameter

Part Number	Color	Length (mm)
126.2036	Magenta	35
126.2041	Magenta	40
126.2046	Magenta	45
126.2051	Magenta	50
126.2056	Magenta	55
126.2061	Magenta	60

SEURAT™ Screwbody: 8.5mm Diameter

Part Number	Color	Length (mm)
126.2241	Light Blue	35
126.2246	Light Blue	40
126.2251	Light Blue	45

